# Cefar TENS

## User Manual

## CONTENTS

<b>1. Introduction</b> .....	4
<b>2. How to use the medical equipment (Intended Use)</b> .....	5
2.1 Fields of Application.....	5
2.2 Therapy objectives.....	5
2.3 Contraindications.....	6
<b>3. Warnings</b> .....	7
<b>4. Overview</b> .....	13
4.1 Device components and accessories.....	13
4.2 Explanation of symbols .....	13
<b>5. Description of the Cefar Tens</b> .....	15
5.1 Description of the device.....	15
5.2 Explanation of symbols.....	16
5.3 Battery Charge.....	17
5.4 Accessories.....	18
<b>6. Step by step use - Patient mode</b> .....	19
6.1 Attach the electrodes.....	19
6.2 Switch on the stimulator.....	19
6.3 Select a preset Program.....	20
6.4 Adjust timer settings if needed .....	20
6.5 Start the stimulation.....	21
6.6 Key Lock.....	21
6.7 Program pause.....	22
6.8 Open lead.....	22
6.9 Resume treatment after open lead.....	23
6.10 Program end.....	23
<b>7. Physician Mode</b> .....	24
7.1 To access Physician Mode.....	24
7.2 Physician Mode home screen .....	24
7.3 Lock Program.....	25
7.4 Custom Programs.....	25
7.4.1 Select custom Programs.....	25
7.4.2 Custom Program first screen.....	25
7.4.3 Summary screen.....	26
7.4.4 Parameters range.....	26
7.5 Compliance.....	27
7.6 CC-CV.....	27

7.7 Accupuncture.....	27
7.7.1 Precautionary measures.....	27
7.7.2 Recommendations.....	28
7.8 Pain score.....	29
<b>8. Programs.....</b>	<b>30</b>
8.1 Program parameters.....	30
8.2 Program description.....	31
<b>9. Troubleshootings.....</b>	<b>33</b>
9.1 Electrode or cable failure.....	33
9.2 Battery level.....	33
9.3 Others.....	34
<b>10. Care, Maintenance, Transport, Environmental Statement.....</b>	<b>36</b>
10.1 Maintenance.....	37
10.2 Transport.....	37
10.2.1 Transport of the CEFAR TENS.....	37
10.3 Environmental Statement, Expected Life.....	38
<b>11. Technical Datas, Standards, Guarantee, Patients.....</b>	<b>39</b>
11.1 Technical Data.....	39
11.1.1 General information.....	39
11.1.2 Stimulation Parameters.....	39
11.1.3 Information on electromagnetic compatibility (EMC).....	40
11.1.4 Environmental conditions.....	40
11.2 Standards.....	41
11.3 Guarantee.....	41
<b>12. Frequently Asked Questions.....</b>	<b>42</b>
<b>13. Contacts.....</b>	<b>43</b>
<b>14. Electrode placement guide.....</b>	<b>44</b>

## 1. INTRODUCTION

As part of Chattanooga, CEFAR stimulators are widely used and endorsed by medical professionals around the world.

TENS is effective and economical. Through clinical research, areas of application for Transcutaneous Electrical Neurostimulation (TENS) and Percutaneous Neurostimulation (PENS) are rapidly expanding. Chattanooga/CEFAR is working actively to further develop the method towards a natural treatment alternative for the patient.

CEFAR TENS is a dual channel nerve stimulator (TENS) designed for use by both professionals and home users. It provides relief of a wide range of pains, like back pain and neuropathic pain. The stimulator features 10 preset programs and 5 custom programs. Program and amplitude can be individually set for each channel, which means that two different programs can be combined simultaneously.

Chattanooga is the world's largest manufacturer of rehabilitation and pain relief equipment. For over seven decades, Chattanooga products have contributed to treating patients in hospitals, clinics and home settings worldwide.

Chattanooga is a brand of DJO Global Inc. DJO Global provides solutions for musculoskeletal health, vascular health and pain management. Our products help prevent injuries or rehabilitate after surgery, injury or from degenerative disease and relieve pain enabling patients to regain or maintain their natural motion. Visit [www.DJOglobal.eu](http://www.DJOglobal.eu)

## 2. HOW TO USE THE MEDICAL EQUIPMENT (INTENDED USE)

### Note

- This manual is considered as an accessory of the therapy unit and therefore it should accompany it at all times.
- The specific instructions provided here are conditions for the intended use and correct operation of the equipment as well as the safety of the patient and the operator using it.

### 2.1 Fields of Application

The CEFAR TENS is a stimulator used by patients at home to ensure electric stimulation treatments in pain management (TENS) and neuro muscular stimulation (EMS/NMES).

The CEFAR TENS is a stimulator used by medical professionals to apply Electro-Acupuncture and neuro muscular stimulation (EMS/NMES) treatments

TENS = Transcutaneous Electrical Nerve Stimulation. TENS is an effective method offering pain relief. TENS uses the central nervous system's own pain relief mechanisms to reduce both acute and chronic pain.

The CEFAR TENS unit is an important supplement to medical and therapeutic treatment for use in hospitals, clinics, general practices and in the patient's home.

### 2.2 Therapy objectives

As a TENS device, indications are for the following conditions:

- Symptomatic relief and management of chronic, intractable pain
- Adjunctive treatment for post-surgical and post-trauma acute pain
- Relief of pain associated with arthritis.

With the Pain Relief programs, impulses are sent through the nerves to inhibit the pain impulses. Pain relief is usually most efficient during stimulation, but the effect can also last after the program has finished. It also increases the circulation. If the pain persists, seek medical advice.

As an NMES device, indications are for the following conditions:

- Retarding or preventing disuse atrophy
- Maintaining or increasing range of motion
- Re-educating muscles
- Relaxation of muscle spasms
- Increasing local blood circulation

The CEFAR TENS is a multifunctional electrotherapy unit that can be used on any occasion when pain relief or muscle relaxation is needed. All the programs can be used as often as required. Guidelines for using the preset treatment programs with indication examples, suggested electrode placements and recommended treatment time are included in this manual.

## 2.3 Contraindications

### **Danger!**

Do NOT use the CEFAR TENS on patients with:

- Implanted electronic devices. Do not use the device if you have a cardiac stimulator, implanted defibrillator or other implanted electronic/electrical device. This could cause an electrical shock, burns, electrical interference or death.
- Epilepsy
- Pregnancy (do not use on abdominal region)
- Serious arterial circulation problems in lower limbs
- Do not use chest stimulation on patients with cardiac arrhythmia or angina pectoris
- Heart disease: If you have suspected or diagnosed cardiopathy you should follow the precautions for use recommended by your doctor because the device may cause lethal rhythm disturbance to the heart in susceptible individuals.

### **Note**

#### **Fixation metalwork and/or prosthesis equipment**

The presence of fixation metalwork and/or prosthesis equipment (metallic equipment in contact with the bone: pins, screws, plates, prostheses, etc) is not a contraindication for TENS and NMES stimulation. The electrical current of the CEFAR TENS is specially designed to have no harmful effect on osteosynthesis equipment.

Electro-stimulation should not be applied if any part of the implanted element going through the patient's skin because this may cause shocks or burns.

### 3. WARNINGS

#### Definitions

It is mandatory to read the safety statements before using the physiotherapy unit. The safety statements are classified as follows:

**Danger!**

This term indicates an imminent hazard. If not avoided, this hazard could result in death or serious injury.

**Warning!**

This term indicates a hazard. If not avoided, this hazard can result in death or serious injury.

**Caution!**

This term indicates a potential hazard. If not avoided, this hazard can result in minor personal injury and/or product/property damage.

#### Safety Information

**Danger!**

Explosion hazard - CEFAR TENS is not designed for use in areas where an explosion hazard may occur. An explosion hazard may result from the use of flammable anaesthetics, oxygen-rich environments, skin cleansing agents and disinfectants.

**Warning!**

Patient hazard

- Only authorised individuals are allowed to repair the CEFAR TENS. Individuals are authorised after receiving training in the operation of the unit and reading this operating into the manual.
- Before using the therapy unit, the operator must ascertain that it is in correct working order and operating condition. The cables and connectors, in particular, must be checked for signs of damage. The device cannot be used if any part is damaged as there is a risk of electrical shock.
- Stop therapy immediately if you have doubts about the device settings and/or the

therapy protocol.

- Patients must be fully conscious while being instructed in the use of the therapy unit and during therapy.
- The choice of the therapy parameters to program and of the therapy protocols to use is restricted to the responsible physician or therapist. It is the physician's or therapist's decision whether or not to use the unit on a specific patient.
- The patients must be familiar with the functions of the CEFAR TENS device allowing them to stop therapy, if needed. Patients unable to operate the emergency stop function e.g. paralytic patients, must never be left unattended during therapy.
- Any accessories used with CEFAR TENS must first be approved by the manufacturer. Using of not specified accessories and detachable parts could be unsafe for the user.
- The utmost caution is advised under the following conditions. Depending on the judgement of the responsible physician, the unit may only be applied under supervision and with the parameters defined by the responsible physician. Otherwise the exercise may be too strenuous for the patients with:
  1. Hypertension (> stage 2), ischemic heart disease and cerebrovascular diseases
  2. Cardiovascular diseases
  3. Pregnancy

Never apply the electrodes:

- Near the head or on the eyes
- On the front and side of the neck (especially carotid sinus) because this could cause severe muscle spasms resulting in closure of the airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure.
- Counter-laterally, i.e. do not use two poles connected to the same channel on opposite sides of the body.
- On or near skin lesions of any kind (swelling, burns, irritation, eczema, cancerous lesion, etc)
- Crossing over the heart because the introduction of electrical current into the chest might cause rhythm disturbance to the patient's heart, which could be lethal.
- If the person is pregnant do not place electrodes directly on the uterus area or connect pairs of electrodes on either side of the abdomen to avoid any risk for the mother and/or the baby. - Keep the device out of the reach of children

### Warning!

- Keep the device out of reach of children
- Never leave the device unattended when it is switched on!
- After use, store the device in a safe place to avoid other people not informed to use the device.
- This device is not a toy but a medical device. Misunderstanding its use can cause damage!

**Warning!**

**Shock hazard** - Strictly observe the following warnings. Failure to do so could endanger the lives of the patient, the user and other persons involved.

- **Before use** allow the CEFAR TENS to reach room temperature. If the unit has been transported at temperatures below 0 °C (32°F), leave it to reach at room temperature for about 2 hours, until any condensation has disappeared.
- **Electrosurgical equipment or defibrillators.** Disconnect the electrodes from the device before using electrosurgical equipment, or a defibrillator, to avoid cutaneous burns from the electrodes and destroying the device.
- **Electronic surveillance equipment.** Do not apply stimulation near electronic surveillance equipment (e.g. cardiac monitors, ECG alarms), as there is a risk they may not work properly whilst the electrical stimulation device is being used.
- **Electromagnetic radiation.** Do not use the stimulator in areas in which unprotected devices are used to emit electromagnetic radiation. Portable communications equipment can interfere with the device.
- **Cancer.** Do not apply stimulation if you have progressive cancer or near any cancerous tumour. The increased metabolism, caused by certain modes of stimulation, is likely to encourage cancer cells to spread
- **Muscle shortening.** During the muscular contraction phase it is recommended to hold the extremities of the stimulated limbs to avoid any shortening of the muscle during contraction, which could cause cramps.
- **Loss of sensation.** Proceed with caution if stimulation is applied to areas of the skin whose level of sensation is lower than normal. Do not apply stimulation to a person who cannot express themselves.
- **Battery leakage.** If there is leak from a component, take steps to ensure the liquid does not come into contact with skin or eyes. Should this occur, wash the affected area with water and consult a doctor.
- **Strangulation.** Do not wind cables around the neck. Tangled cables can cause strangulation.
- **Post-surgery.** Proceed with caution after recent surgery when stimulation may disrupt the patient's healing process.
- **Accessibility of the power adaptor.** The plug socket must be close to the power adaptor and be easily accessible.
- The CEFAR TENS must only be operated in dry rooms.
- Do not use the CEFAR TENS in water or in a humid atmosphere (sauna, bath, shower etc.), etc.).
- When connecting the unit to other equipment or when creating a medical system, check that the sum of leakage currents will not cause any hazard. Please contact DJO GLOBAL if you have questions regarding this matter.
- No modification of this equipment is allowed.
- Do not open the product and its accessories as there is risk of electrocution.
- Before cleaning and service interventions, **turn off the device.**

- **Liquids and foreign material (such as dust, metal etc) must not be allowed** to enter the device. If such material has entered into the units, it must be immediately checked by a service technician, before it can be reused.
- **Electricity supply.** Never connect the stimulation cables to an external power supply, as there is a risk of electrocution.
- Do not apply stimulation near the area of an implant, such as cochlear implants, pacemakers, skeletal anchorage or electric implants. This could cause an electrical shock, burns, electrical interference or death.
- Never use the CEFAR TENS or the AC adaptor if it is damaged or open. There is a risk of electric shock.
- Disconnect the AC adaptor immediately if there is abnormal heating or smell, or if smoke comes from the AC adaptor or the device.

### Warning!

Equipment malfunction - this warnings can prevent equipment malfunctions that result in patient hazards

- Keep the unit away from such equipment and verify its performance before use.
- Do not use the CEFAR TENS within one metre of short wave or microwave devices as this could alter the currents generated by the stimulator. If you are in any doubt as to the use of the stimulator in close proximity to another medical device, seek advice from the manufacturer of the latter or from your doctor.
- Exercise caution when using electrotherapy while the patient is connected to monitoring equipment with electrodes attached to the body. Stimulation could disrupt the signals sent to the monitoring equipment.
- Refer repair and maintenance to authorized persons. Persons are authorized after training by a specialist trained and commissioned by the manufacturer.
- Inspect the CEFAR TENS and it's accessories for damage and loose connections at least once a year. Damaged and worn parts must be immediately replaced with original spare parts by authorized staff.
- Do not use the device if you are connected to a high-frequency surgical instrument as this could cause skin irritation or burns under the electrodes.

### Caution!

Patient hazard - these cautions need to be observed to avoid the risk of electrical shock or other negative effects to the patient.

- Do not apply stimulation close to metal. Remove jewellery, piercings, belt buckles or any other metallic product or device in the area of stimulation
- Be careful if the patient has sensitivity problems or is not able to communicate that he

- or she feels discomfort, however light.
- Never begin an initial stimulation session on a person who is standing. The first five minutes of stimulation must always be performed on a person who is sitting or lying down. In rare instances, people of a nervous disposition may experience a vasovagal reaction. This is of psychological origin and is connected with a fear of the muscle stimulation as well as surprise at seeing one of their muscles contract without having intentionally contracted it themselves. A vasovagal reaction causes heart to slow down and blood pressure to drop, which produces a feeling weakness and a tendency towards fainting. If this does occur, all that is required is to stop the stimulation and for the person to lie down with the legs raised until the feeling of weakness disappears (5 to 10 minutes)
  - Do not disconnect any channel during the stimulation.
  - Do not use the stimulator while driving or operating machinery or during any activity in which electrical stimulation or involuntary muscle contraction can put the patient at risk of injury.
  - Do not apply stimulation during sleep.
  - Do not use the stimulator at altitudes of over 3,000 meters.
  - When using constant current mode, always turn off the stimulator before moving or removing any electrodes during a session, to avoid electrical shock to the patient.
  - Do not try to place electrodes on a body part not directly visible without assistance.
  - Attach the electrodes in such a way that their entire surface is in contact with the skin.
  - For obvious reasons of hygiene, each patient must have their own set of electrodes. Do not use the same electrodes on different patients.
  - Some patients with very sensitive skin may experience redness under the electrodes after a session. Generally, this redness is completely harmless and usually disappears after 10 to 20 minutes. Never start another stimulation session in the same area, however, if the redness is still visible.
  - When using the device to customize programs, take special care that the parameters customized and applied by you to the patient are as you wanted them to be.

### Caution!

- Do not expose the CEFAR TENS to direct sunlight, because some of the components may reach unacceptably high temperatures.
- Make sure that pets or vermin do not contaminate the physiotherapy unit and keep them away from it. Also, keep the unit clean and protect it from dust and lint. The safety rules and regulations set forth apply in any case.

**Caution!**

- It is recommended to use the transport bag that comes with the unit, for transport of the device, and to use a proper transport box to ship it.
- Always use the AC adaptor (power supply) provided by the manufacturer to recharge the unit. Never use any other AC Adapter. This could damage the device and generate electrical hazards.
- Only use electrodes supplied by the manufacturer. Other electrodes may have electrical properties that are unsuitable for or may damage the CEFAR TENS.
- Do not place the electrodes in water.
- Do not apply solvents of any kind to the electrodes.
- **Instructions for electrodes.** See the usage and storage instructions displayed on the bag of electrodes.
- **Hot casing or battery.** Under extreme use conditions, some parts of the casing might reach up to 118 °F (48 °C). Use caution when manipulating the rechargeable battery right after device use or when holding the device. There is no particular health risk associated with this temperature besides your comfort
- **Size of electrodes.** Do not use electrodes with an active area of less than 10 cm<sup>2</sup> due to the risk of associated burning. Proceed systematically with caution when the density of the current is over 2 mA/cm<sup>2</sup>.
- If the device is not used for an extended period, we recommend that you charge your device every 5 months. You should store the device in a dry and cool environment.

**Note**

- For best results, wash and clean the skin of any oil and dry it before attaching the electrodes.
- Never use a set of adhesive electrodes for more than 40 sessions as the quality of the contact between the electrode and the skin, which is essential for the patient's comfort and the effectiveness of the stimulation, gradually deteriorates.
- For information on use and storage please consult the instructions found on the electrodes packaging.

**Note****Biocompatibility**

Those parts of the CEFAR TENS unit that come into contact with the patient when used as intended, are designed to fulfil the biocompatibility requirements of the applicable standards.

## 4. OVERVIEW

### 4.1 Device components and accessories

**Model:** CEFAR TENS

**Part number:** 2053/2061/2075/2094/2106

Your kit contains (included in delivery):

Quantity	Description	Part number
1	CEFAR TENS device	1981
1	Fast charger	2001 (EU), 2133 (UK), 2145 (US)
1	Set of 2 pin cables 2 mm	1990
1	User manual and practical guide	13-00187
1	Transportation pouch	2049
2	Necklace and clip	2023
2	Beltclip	2037

Electrodes are considered as option depending on markets.

### 4.2 Explanation of symbols (connections and nameplates)




	Follow the instructions of the user manual
	Caution! Observe warnings set forth in operation manual!
	The CEFAR TENS is a class IIa device with internal electric power and type BF applied parts.
	The power switch On/OFF button is a multi-function button.
	The name and address next to this factory symbol is the manufacturer.

	The number next to this symbol is the article reference number
CE-2797	Device complies with Council Directive 93/42/EEC as amended, concerning medical devices.
	The number next to this symbol is the serial number
	WEEE Mark (European Directive 2002/96/EC). Indicates separate treatment from general waste at end of life.
	Keep dry
IP22	IP classification indicates the degree of protection and thus defines its suitability for use under various ambient conditions.
	Non-ionising radiation
	Not made with natural latex rubber
	Keep away from direct sunlight
	Alternating current input on AC power supply
	Direct current output from power supply
	Protection class II equipment. The AC Power Supply device has double insulation.
	Indoor use only

## 5. DESCRIPTION OF THE CEFAR TENS






### 5.1 Description of the device

- A ON/OFF
- B Program selection - Pause during stimulation
- C Increase/decrease amplitude increase/decrease timer menu navigation
- D Timer
- E Confirm - Resume stimulation when paused
- F Cables sockets
- G Battery charger socket
- H Belt clip holders
- I Belt clip (can be reversed 180°)

## 5.2 Explanation of symbols

- 1 Battery Status
- 2 Keyboard lock symbol
- 3 Program number (left and right channel)
- 4 Amplitude bargraph (left and right channels)
- 5 Amplitude level in milli Ampere (mA) or Voltage (V)
- 6 Remaining time
- 7 Mixed freq - Muscle stimulation icon
- 8 Pause
- 9 Program end
- 10 Open circuit (disconnect electrode or cable)
- 11 Program Lock

### 5.3 Battery charge

The CEFAR TENS is powered by a rechargeable battery.

- A - To recharge the CEFAR TENS first disconnect the electrode cables from the device, then plug the charger to a wall socket then connect the stimulator to the charger.
- B - The battery level of charge is indicated on the treatment screen
- C - When charging, the battery charge symbol is animated with a battery filling itself from empty to full.
- D - When battery is fully charge at the end of charging process device will show 100%.

Fig A


Fig B


Fig C - Low charge


Fig D - Half charged


Fig E - Full charge

### 5.4 Accessories

**Electrodes:** The CEFAR TENS should be used with recommended DJO electrodes. The electrodes will eventually wear out and need to be replaced. It is recommended to replace after approximately 20-40 times of usage.

**Note:** Stimulation comfort and treatment effectiveness are highly influenced by the quality of the electrodes.

**Neck Strap and Clip:** The CEFAR TENS comes with a belt clip and a neck strap/clip allowing you to wear the stimulator around your neck or on your belt and have your hands free during treatment.

To use the neck strap, fix the clip associated to the hole at the top of the unit. (Fig A & Fig B)

The Belt clip can be rotated in order to have the cable outputs on the bottom, or on the top depending on the body area being stimulated. To release the belt clip, push it downwards (Fig C & Fig D).

Fig A

Fig B

Fig C

Fig D

## 6. STEP BY STEP USE - PATIENT MODE

You must read the contra-indications and safety measures described at the beginning of this manual (Chapter 3 "Safety Information") prior using your stimulator.

### 6.1 Attach the electrodes

- Connect the electrodes to the cables.
- Attach the electrodes to your body.
- Plug the stimulator cables into the two sockets on the base of the device. The two cables can be connected simultaneously. The cables colours are different to simplify use and facilitate identification of the two channels. Make sure to have the line of the connector upwards to insert the connector into the right direction.
- The kit includes two 2mm pin-connector cables.

### 6.2 Switch on the stimulator

Press the "On/Off" button.

**Note:** last program configuration will always be displayed by default.

### 6.3 Select a preset Program

Press the Program "P" button for the left or right channel you want to adjust. The program number is blinking on the screen. Push on up/down arrows of the same channel to increase or decrease the program number and confirm with "OK". If you want to use both channels, then select a program for the other channel. Any combination of programs is possible except for Program 7, which must be selected for both channels.

**Note:** When selecting a program the amplitude must be set to 00.0 mA or 0.0V for the channel concerned. For further information on the programs, see section PRESET PROGRAMS.

### 6.4 Adjust timer settings if needed

Timer preset value can be either unlimited (---) or with a limited time. See program table summary at the end of this manual. Press once the "timer" button, Channel 1 time value is blinking. To adjust time press on adjust value buttons "up or down" and press "OK" to confirm

Press twice the timer button to select Channel 2 timer adjustment.

## 6.5 Start the stimulation

Press the "INCREASE" button for each channel until you reach a comfortable level of stimulation. Press and hold the button to increase the amplitude gradually.

**Note:** Always increase the amplitude cautiously!

**Note:** Program 5 (TENS HAN): Set the amplitude for conventional stimulation (80 Hz) when the upper part of the high/low frequency symbol is flashing. When the stimulator switches over to burst after 3 seconds of stimulation (2 Hz) and the bottom part of the symbol is flashing, set the amplitude for this frequency. The frequencies switch every 3 seconds.

## 6.6 Key Lock

An automatically activated key lock prevents unintentional changes during treatment. The key lock is activated if no button is pressed for 10 seconds while stimulation is ongoing. Press any of the "DECREASE" buttons to deactivate the key lock.

## 6.7 Program pause

To pause a program press "P", both Channels 1 and 2 are paused, Amplitude goes down automatically to 00,0. To resume program press "OK". Amplitude increases automatically and progressively to set value.

## 6.8 Open lead

The broken circuit display indicates on which channel that resistance is too high, or that a cable is broken or detached. A high resistance can be caused by a bad connection between the electrodes and your skin, or that the electrodes need to be replaced. The stimulation is stopped, while trouble shooting is fixed.

**Note:** Broken circuit is only indicated when running over 9.5mA. Arrow is pointing to the channel with the open lead.

### 6.9 Resume treatment after open lead

Press "OK" after checking electrodes and cables connections. Treatment screen is displayed. Amplitudes have to be readjusted manually. A cable breakage can be checked by pressing the cable's pins against one another while increasing the amplitude for the corresponding channel to 11 mA. If the amplitude now drops to 0.0 mA and starts flashing, the cable needs to be replaced.

**Note:** Never increase the amplitude above 20 mA when you check for cable breaks, since this can damage the stimulator.

### 6.10 Program end

When the program is finished, the last used programs are selected for the next time the stimulator is turned on.

The stimulator turns off after 5 minutes of inactivity to spare battery.

## 7. PHYSICIAN MODE

You must read the contra-indications and safety measures described at the beginning of this manual (Chapter 3 "Safety Information") prior using your stimulator.

Physician mode can only be used by trained medical professionals. Other people should not use this mode.

### 7.1 To access physician mode

When device is switched OFF, press simultaneously "Ch2+" and "P2" buttons and maintain them pressed while pushing on the "ON/OFF" button.

### 7.2 Physician mode home screen

Navigation using "up and down", press "OK" to confirm

**Lock Program:** The CEFAR TENS can be locked to prevent accidental changing of programs

**Custom Program:** With the CEFAR TENS it is possible to create and store five custom programs for individual needs

**Access Compliance:** Compliance gives you the possibility to monitor the use of the stimulator

**Reset Compliance:** Reset compliance data

**Set CC-CV:** CEFAR TENS can operate on Constant Voltage or Constant current mode

**Select Acupuncture:** The CEFAR TENS offers a dedicated electrical acupuncture mode: Acupuncture mode can only be used by trained medical professionals.

**Enter Pain score:** Patient pain can be assessed with numerical scale

**Reset pain score:** Patient pain score can be reset

### 7.3 Lock Program

It is possible to lock a program to avoid using an incorrect program. After selecting a program in patient mode, use up and down to select "Lock/Unlock Program", "OK" to confirm and go back to Physician home screen.

### 7.4 Custom Programs


#### 7.4.1 Select Custom Programs

P11;P12;P13;P14;P15 can be customized for individual needs.

Use "up and down" to select program, "OK" to confirm.

#### 7.4.2 Custom Program first screen

Stimulation mode can be set as Conventional, Burst, Pulse Width Modulation (PWM), Frequency Modulation, Mixed Frequency. Press "up and down" to select appropriate mode in the parameters range. Press "OK" to confirm.


### 7.4.3 Summary screen

Summary screen is displayed after last parameter is adjusted and "OK" button is pressed once. Use "up and down" to scroll through the different parameters for visualization only.


<b>P 12</b>	
<b>Conventional</b>	
Frequency	<b>80</b>
Pulse duration	<b>180</b>

### 7.4.4 Parameters range

### 7.5 Compliance

Device displays session time as time of use since last reset. Total time of device use can't be reset. Session time can be reset.



### 7.6 CC-CV

CEFAR TENS can be set on Constant Current or Constant Voltage mode. Press "up and down" to select appropriate mode. "OK" to confirm and back in Physician home screen.

**CC - Constant current:** This mode guarantees optimal stimulation without variation of amplitude by automatically adapting the voltage when resistance varies during the stimulation. The resistance may vary depending on tissue impedance, or electrodes contact quality. The CC mode is recommended by default.

**CV - Constant Voltage available only on 1 channel:** This mode adapts automatically the amplitude when the resistance varies, allowing more comfortable stimulation if the patient loses electrode surface contact or takes electrodes the electrodes in their fingers while stimulation is on. CV mode is recommended to be used when moving electrodes from place to place or when the patient does not feel comfortable managing the device.

### 7.7 Acupuncture

Acupuncture mode can only be used by trained medical professionals. Other people should not use this mode. Medical professional must constantly be with the patient while giving the treatment.

#### 7.7.1 Precautionary measures

The Acupuncture Electrical stimulation can only be applied by authorized medical professionals using dedicated electrical acupuncture needles. When Acupuncture mode is selected the "Acupuncture" icon is displayed on the treatment screen. Using Acupuncture mode the amplitude can be set from 0,1mA steps. Be careful when stimulating patients with a nickel allergy.

Note that significant effects can be obtained with electrical acupuncture-overstimulation may entail worsening of the patient pain and general condition. When stimulating at less than 10 mA in acupuncture mode, the stimulator doesn't register that the circuit is broken – which means that the current is not disconnected even if one of the clamps comes loose. Always make sure that the amplitude is turned off before connecting the clamps.

## 7.7.2 Recommendations

Select Acupuncture ON/OFF in the physician mode.

Acupuncture
<b>ON</b>
OFF

Do not place needles for electro-acupuncture from the same pair directly over the heart.

Disinfect and dry the skin before penetrating with the acupuncture needles. Connect the clamps to the needles. Place the clamps as close to the skin as possible to prevent "sagging" needles.

The distance between the two needles in a pair should not be less than 3 cm (1- 1/4") or greater than 30 cm (12"). The patient must feel the stimulation clearly in the treatment area.

Both high-frequency and low-frequency stimulation should be felt clearly, but with absolutely no pain. With high-frequency stimulation, the patient should clearly feel paresthesia, and with low-frequency stimulation, visible muscle contractions should occur. This can sometimes be difficult to attain without painful stimulation (for example in the facial area). If so it is important to clearly see muscular fasciculation around the needle.

### High-frequency stimulation:

Gives segmental, fast, brief effect. Used primarily with local, segmental needles. Usually acute and sub-acute pain conditions are treated with high-frequency stimulation for rapid pain alleviation.

### Low-frequency stimulation:

Gives segmental and general effect, slower. Used with local and distal extra segmental needles. Used primarily to amplify the general effect of manual acupuncture.

### Mixed frequency stimulation:

Gives a combined effect of low-frequency and high-frequency stimulation.

Note when treating with mixed-frequency stimulation: You must set amplitude for each frequency.

**Stimulation with short pulse duration:**

The pulse duration is a way to regulate the amount of energy. A small pulse duration is suitable for treating sensitive, nerve-rich areas like the face.

**7.8 Pain score**

- Pain score can be assessed and stored in a numerical scale from "0" absence of pain to "10" maximum pain.
- Use "**Up and Down**" to set pain score
- Press "**OK**" to confirm
- Pain score can be reset

Pain score (0-10)  <b>5</b>	Pain score <b>5</b> <b>Reset</b> Cancel
-----------------------------------	---

## 8. PROGRAMS

### 8.1 Program parameters

Program Number	Name	Pulse Width 1	Pulse Width 2	Frequency 1	Frequency 2
P1	Gate Control HF 80 Hz	180 $\mu$ s		80 Hz	
P2	Burst LF 2 Hz	180 $\mu$ s		2 Hz	80 Hz
P3	Gate Control HF 80Hz Mod	70 $\mu$ s	180 $\mu$ s	80 Hz	
P4	Gate Control HF 80Hz Sensitive	60 $\mu$ s		80 Hz	
P5	Mix 80 Hz/2 Hz Han	180 $\mu$ s		80 Hz	2 Hz
P6	Neuromodulation LF 10 Hz	180 $\mu$ s		10 Hz	
P7	HF Flow Tens	70 $\mu$ s	180 $\mu$ s	80 Hz	
P8	Endorphinic LF 5 Hz	180 $\mu$ s		5 Hz	
P9	Gate Control HF 100 Hz	180 $\mu$ s		100 Hz	
P10	NMES Muscle Stim	250 $\mu$ s		50 Hz	0 Hz
P11	Custom 1				
P12	Custom 2				
P13	Custom 3				
P14	Custom 4				
P15	Custom 5				

## 8.2 Program description and recommendations

P1	<p>Conventional TENS (high-frequency stimulation) is a good choice for both acute and long-term pain, both neurogenic and nociceptive. Conventional TENS is based on the Gate Control theory, which states that electric stimulation of A-beta fibers inhibits impulse transfer in the pain pathways. As a rule, the electrodes should be placed on or near the painful area, or over an area segmentally related to the painful area. Adjust the amplitude so that the stimulation gives strong, but pleasant paraesthesia – tingling. When using high-frequency stimulation, it's important to ensure that the patient has normal sensory of touch in the area where the electrodes are placed.</p>
P2	<p>Burst TENS (low-frequency stimulation) is usually most effective for radiating (projected) pain in the arms and legs (rhizopathy), for conditions with reduced or changed sensory of touch, for deep muscular pain, or when the post-treatment effect of conventional TENS is too short. Burst TENS treatment alleviates pain by stimulating muscles to release the body's own morphine-like substances, endorphins. Place the electrodes on a muscle in the painful area so that a visible contraction occurs, or on acupuncture points in the painful area. The stimulation should feel pleasant and give visible muscle contractions. Remember that the patient often feels the stimulation clearly before contractions become visible.</p>
P3	<p>Modulated pulse duration stimulation is a type of high-frequency stimulation where the pulse duration varies continuously. This can cause an undulating sensation, which may be more pleasant than a constant pulse duration. Use Program 3 for pain alleviation and a massage effect on muscles like the trapezius.</p>
P4	<p>Program 4 has a shorter pulse duration than the other programs. A short pulse duration is suitable for treating sensitive, nerve-rich areas like the face and high up on the neck. With a shorter pulse duration the amplitude can be increased, which makes it easier to find the right level without the patient experiencing pain.</p>
P5	<p>Mixed frequency stimulation is also called Han stimulation (3 sec at 2 Hz and 3 sec at 80 Hz). Stimulation frequencies switch every three seconds, giving a combination of high- and low-frequency stimulation, which can offer a more effective treatment. Place the electrodes as you would for low-frequency stimulation – on a muscle in the painful area.</p>
P6	<p>Program 6 is specifically designed for treating acupuncture points on a neuromodulation mode. For nausea, place the black electrode on acupuncture point PC6 (inside of the wrist) and the red one on acupuncture point TE5 (outside of the wrist). For pudendal pain as for incontinence treatment, place the electrode on SPE nerve on medial part of the ankle.</p>


P7	Always use 2 channels and 4 electrodes with Program 7. A modulated pulse duration stimulation is obtained, but since the channels are alternately active, a massaging, pumping effect is experienced. Place the electrodes as for high-frequency stimulation when a pain relieving and massaging effect is desired.
P8	An increase in the tension of the contracted muscle fibres and the crushing of the capillary network resulting from this causes a decrease in the blood flow and a gradual accumulation of acid metabolites and free radicals. Without treatment, there is a risk that the contracture will become chronic and genuine atrophy of the capillary network may gradually occur. In addition to the general effect of increasing endorphin production in the hypothalamus, which elevates the pain perception threshold, there is a very significant localised effect. The 5 muscle twitches induced every second by stimulation produce very significant hyperaemia, which drains the acid metabolites and free radicals that had accumulated in the chronically contracted muscle areas.
P9	The principle is to cause a significant influx of tactile sensitivity in order to restrict the entry of pain impulses upon their return to the posterior horn of the spinal cord. We must therefore stimulate the sensitivity fibres on the skin of the painful area. To do this, it is necessary to use a frequency that is the same as the operational frequencies for the tactile sensitivity nerve fibres that could attend 100 Hz
P10	A muscle that is normally innervated, after a period of immobilisation or diminished movement due to pain or restricted range of motion, rapidly decreases in volume. This decrease depends on the degree and duration of the functional deficit. Slow fibres (type I) in particular are affected by disuse atrophy.
P11-P15	In the Programming mode you could create a custom program and select your own parameters

## 9. TROUBLESHOOTING

### 9.1 Electrode or cable failure

The CEFAR TENS will display a pair of electrodes and an arrow pointing to the channel where the problem was detected.

The electrode error message can mean that:

- There are no electrodes connected to this channel.
- The electrodes are old, worn, and/or the contact is poor: try using new electrodes.
- The electrode cable is faulty: try connecting it to another channel. If the problem persists, replace the cable.

If for any reason your device seems to present a malfunction other than those presented above, please contact DJO Global customer service.

Problem	Possible Cause	Solution
<b>Electrode failure</b>	Poor electrode connection to the cable	Check that the electrodes are properly connected to the cable.
	Poor electrode connection to the skin	Check if the electrodes are outdated, worn and/or the contact is poor: try with new electrodes

### 9.2 Battery Level

Problem	Possible Cause	Solution
<b>Device battery low</b>	During the stimulation the device may turn out due to discharged battery.	Stop the stimulation and recharge the device.

## 9.3 Others

Problem	Possible Cause	Solution
Display does not come on	Low battery	Charge battery
Weak stimulation with loaded batteries	Electrodes dried out, lost their adhesive power and have no adequate connection to the skin.	Replace electrode
	Acupuncture mode is wrongly selected	Unselect acupuncture mode
	Electrode placement	Make the electrodes at least 5 cm (2 inches) apart
Stimulation stops with loaded batteries	Poor electrode contact	Reapply electrodes, secure firmly. Electrodes must be a minimum of 5 cm apart.
	Damaged or worn electrodes	Replace
	Cabel or contacts damaged	Replace
Stimulation weakens within minutes of starting treatment with loaded batteries	This is a normal body adaptive process	Increase the amplitude (intensity) if required.
Stimulation is uncomfortable	Amplitude (intensity) is too high	Decrease amplitude (intensity).
	Electrodes are too close together	Reposition the electrodes. Electrodes must be a minimum of 5 cm (2 inches) apart.
	Damaged or worn electrodes	Replace
	Ensure proper program is being used	A. Refer to the programs description B. Contact clinician if discomfort persists

<b>Stimulation is ineffective</b>	Improper electrode placement	Reposition electrodes. Electrodes must be a minimum of 5 cm (2 inches) apart.
	Unknown	Contact clinician.
<b>Stimulation only felt on one electrode</b>	Improper electrode placement	A. Reposition electrode. Electrodes must be a minimum of 5 cm (2 inches) apart. B. Replace electrodes.
<b>Stimulation on one channel (side) only</b>	Electrodes A. Worn or damaged B. Improper placement	A. Replace B. Reposition electrode. Electrodes must be a minimum of 5 cm (2 inches) apart.
<b>Intermittent Output</b>	Intermittent Program in use	Some programs will seem intermittent. This is expected. Refer to section 8.1 for a description of the programs.
<b>Stimulation is not producing the usual sensation</b>	Settings and Electrodes positioning	A. Check that all the settings are correct and ensure the electrodes are positioned properly B. Change the positioning of the electrodes slightly.

When an error is detected by the device (from anywhere in the device), the stimulation shall stop on all channels and a specific screen will be displayed. If this screen is displayed, record the Error number and call the DJO Customer Service.

## 10. CARE, MAINTENANCE, TRANSPORT, ENVIRONMENTAL STATEMENT

### Warning!

Shock hazard - Remove the power cord of the device from the wall outlet before cleaning.

Shock hazard, equipment damage -

- Liquids must not enter the device and its components. If liquids have entered into the components, the device must be immediately checked by a service technician, before it can be reused.
- Never dismantle the device or the AC adapter as they contain high-voltage parts with a risk of electric shock.

- All parts of the device can be disinfected weekly by wiping down with a disinfectant. Thus, it complies with the special hygiene standards for medical technical equipment.
- All components can be cleaned with common disinfectants and mild household detergents.
- Only use a soft cloth and an alcohol-based, solvent-free cleaning product, to wipe the therapy unit down.
- Allow the device to completely dry before use.
- Please refer to environmental part (11.1.4) for storage and transport conditions.

### Warning!

Patient hazard - patient contamination

- Before using the unit on another patient, clean and disinfect it according to the instructions in this section.

### Caution!

Equipment damage -

- The plastic material used is not resistant to mineral acids, formic acid, phenols, cresols, oxidants and strong organic or inorganic acids with a pH value below 4.
- Use only clear disinfectants to prevent discolouration of the device.
- Do not expose the therapy unit to strong ultraviolet radiation (sunlight) and fire.
- Do not sterilize the stimulator.
- Do not immerse in liquids.

## 10.1 Maintenance

Your device does not require calibration or frequently safety testings. Each stimulator is tested prior to distribution. Its characteristics do not vary under normal conditions. If your stimulator contains parts that seem worn or defective, please stop using it and contact the customer service centre that has been stipulated and authorised by the manufacturer regarding an upgrade.

There are no user serviceable parts inside the device. Please contact DJO if you need assistance setting up, using or maintaining the CEFAR TENS or to report any unexpected operations or events. Please also contact DJO for clinical questions or if any component of the CEFAR TENS is not functioning properly or requires servicing. You can find the different contact possibilities under [www.DJOGlobal.eu](http://www.DJOGlobal.eu)

### **Warning!**

Shock hazard, Equipment damage -

Do not attempt to repair the stimulator or any of its accessories. Never dismantle the device because of risk of electric shock. DJO declines all responsibilities for any damages or consequences resulting from unauthorised attempts to open, modify, or repair the stimulator. This may only be done by persons or repair services authorised by the manufacturer.

## 10.2 Transport

### 10.2.1 Transport of the CEFAR TENS

1. Prepare the device and its accessories for shipping within the original CEFAR TENS pouch
2. Turn the device and it's accessories off.
3. Disconnect and dismount the device and it's accessories by following the guidelines
4. Place the accessories within the box
5. Store the user manual in the pocket of the transport bag.

### **Caution!**

Equipment damage -

Only use the original transport bag for carrying the device around.

**Caution!**

Equipment damage -

Only use the original shipping box for shipping the device. DJO cannot be held liable for transport damage if the device is not packed in its original shipping box.

### 10.3 Environmental statement, expected life

The device is electronic equipment and may include substances that can damage the environment. It must not be disposed of with unsorted household or municipal waste. It requires separate disposal at a suitable collection point for recycling of electronic equipment. By doing so, you will be contributing to the safeguarding of natural resources and health. Please contact DJO for information about the possible recycling of the product.

When the electrodes no longer stick well to your skin, dispose of them in a receptacle out of reach of children and pets.

The product as well as the parts and accessories supplied with it are designed for a minimum service life of 5 years of normal usage.

## 11. TECHNICAL DATAS, STANDARDS, GUARANTEE, PATIENTS

### 11.1. Technical data

#### 11.1.1 General information

**Model:** CEFAR TENS

**Part number:** 1981

**MDD:** Class IIa

**IP class:** IP22

**Applied part:** Type BF

**Battery:** Lithium Polymer (Li-Po) rechargeable battery (3.7V 1500 mAh).

**Battery chargers:** Only battery chargers provided by DJO and with P/N 2001 (EU), 2133 (UK), 2145 (US) can be used to recharge the batteries supplied with CEFAR TENS stimulators.

**Weight:** 156g

**Dimensions:** 132 x 60 x 24 mm

#### 11.1.2 Stimulation parameters

All electrical specifications are given for an impedance of 1500 ohms per channel (60mA) and for 1000 Ohms per channel for 99.5 mA.

**Outputs:** Two independent and individually adjustable channels that are electrically isolated from each other.

**Pulse shape:** Constant asymmetric rectangular current with pulse compensation for TENS programs to eliminate any direct current component to prevent residual polarisation at skin level. Symmetric for P10 NMES program.

**Maximum pulse intensity:** 99.5 mA in Constant Current mode or 99 V in Constant Voltage mode for neurostimulation. 30 mA for electroacupuncture.

**Pulse intensity increments:** Manual adjustment of stimulation intensity from 0 to 99 mA in minimum increments of 0.5 mA for Constant Current Mode, 0 to 99 V in 1 V increments for Constant Voltage mode for electrostimulation. 0,1 mA – 30mA steps 0,1 mA for electroacupuncture (5 V for first step).

**Pulse width:** 60 to 400  $\mu$ s.

**Frequency:** 1-120 Hz

**Maximum electrical charge per pulse:** 80 micro coulombs ( $2 \times 40 \mu$ C, compensated)

### 11.1.3 Information on electromagnetic compatibility (EMC)

The CEFAR TENS is designed to be used in typical environments that have been approved in accordance with the EMC safety standard EN 60601-1-2.

This device complies with the CISPR standard, indicating that radio frequency (RF) emissions are not likely to cause interference with electronic equipment installed nearby (radios, computers, telephones, etc.).

The CEFAR TENS is designed to withstand foreseeable disturbances from electrostatic discharge, magnetic fields from the mains power supply or RF transmitters.

Nevertheless, it is not possible to ensure that the stimulator will not be affected by mobile powerful RF (radio frequency) fields from other sources.

For more detailed information concerning electromagnetic emissions and immunity, refer to the EMC tables.

### 11.1.4 Environmental conditions

#### Storage and Transport Conditions

The device must be stored and transported in accordance with the following conditions:

**Temperature:** -20°C to 70°C

**Maximum relative humidity:** 75%

**Atmospheric pressure:** 700 hPa to 1,060 hPa

#### Conditions of use

**Temperature:** 10°C to 40°C

**Maximum relative humidity:** 30% to 75%

**Atmospheric pressure:** 700 hPa to 1,060 hPa

## 11.2 Standards

The CEFAR TENS has been designed, manufactured, and distributed in compliance with the requirements of European Directive 93/42/EEC, as amended, on medical devices.

Current international standards require that a warning be given concerning the application of \ electrodes to the thorax (increased risk of cardiac fibrillation).

The CEFAR TENS also complies with Directive 2012/19/EU on waste electrical and electronic equipment (WEEE).

## 11.3 Guarantee

This guarantee is valid only if it is accompanied by proof of purchase.

Your statutory rights are not affected by this guarantee.

Your CEFAR TENS stimulator is guaranteed for a period of 3 years from the date of purchase. The guarantee covers the device and the AC adaptor (hardware and labour), but not the cables, battery and electrodes.

All defects resulting from poor quality material or workmanship are covered. This guarantee does not cover damage resulting from impact, accidents, misuse, inadequate protection against moisture, immersion in water or repairs made by unauthorised personnel.

## 12. FREQUENTLY ASKED QUESTIONS

### **FOR HOW LONG TIME CAN I STIMULATE?**

Follow the recommendations given by your health care provider. Generally high-frequency TENS (80 Hz) can be used without an upper time limit, but at least 30 minutes at each occasion. Low-frequency TENS (2 Hz) can cause sore muscles but normally 20-45 minutes 3 times a day is recommended.

### **CAN I USE TENS ON ALL KINDS OF PAIN?**

TENS can be used for both acute and long time pain. General instructions can be read in chapter PROGRAM DESCRIPTION. Contact your health care provider if the cause for and area of treatment is changed.

### **WHICH DISTANCE SHOULD I HAVE BETWEEN THE ELECTRODES?**

It is recommended to have a distance of 5 to 30 cm between the electrodes.

### **FOR HOW LONG LAST THE ELECTRODES?**

The self-adhesive electrodes last for approximately 20 to 40 occasions. The durability depends on how good the care and maintenance instructions are followed. Read the care and maintenance instructions delivered with the electrodes.

### **ARE THERE ANY PROOFS OF TENS WORKING AS AN ALLEVIATION OF PAIN?**

TENS is an established and well-documented method that has been used by the medical service for many years.


## 13. CONTACTS

We would be happy to answer any questions you may have about our products and services. Please contact your local dealer, or your corresponding DJO site. DJO sites are listed on the backside of the cover.

For technical service from DJO, please contact:

[internationalproductsupport@djoglobal.com](mailto:internationalproductsupport@djoglobal.com)

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## 15. EMC-TABLES

**Warning!**

Use of the Cefar TENS adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, the Cefar TENS and the other equipment should be observed to verify that they are operating normally.

**Warning!**

Use of accessories and cables other than those specified or provided by the DJO could result in increased electromagnetic emissions or decreased electromagnetic immunity of the Cefar TENS and result in improper operation.

**Warning!**

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Cefar TENS, including cables specified by the DJO. Otherwise, degradation of the performance of the Cefar TENS could result.

Declaration - electromagnetic emission

Cefar TENS is intended for use in the electromagnetic environment specified below. The customer or the user of Cefar TENS should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	Cefar TENS uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	Cefar TENS is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Declaration - electromagnetic immunity

Cefar TENS is intended for use in the electromagnetic environment specified below. The customer or the user of Cefar TENS should assure that it is used in such an environment.

Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.

<p><b>Electrical fast transient/burst</b> IEC 61000-4-4</p>	<p>± 2 kV for power supply lines ± 1 kV for input/output lines</p>	<p>± 2 kV for power supply lines ± 1 kV for input/output lines</p>	<p>Mains power quality should be that of a typical commercial or hospital environment.</p>
<p><b>Surge</b> IEC 61000-4-5</p>	<p>± 0.5kV, ± 1 kV line(s) to lines ± 0.5kV, ± 1 kV, ± 2 kV line(s) to earth</p>	<p>± 0.5kV, ± 1 kV line(s) to lines ± 0.5kV, ± 1 kV, ± 2 kV line(s) to earth</p>	<p>Mains power quality should be that of a typical commercial or hospital environment.</p>
<p><b>Voltage dips, short interruptions and voltage variations on power supply input lines</b> IEC 61000-4-11</p>	<p>0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°and 315°  0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°  0 % UT; 250/300 cycles</p>	<p>0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°and 315°  0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°  0 % UT; 250/300 cycles</p>	<p>Mains power quality should be that of a typical commercial or hospital environment. If the user of the Cefar TENS requires continued operation during power mains interruptions, it is recommended that the Cefar TENS be powered from an uninterruptible power supply or a battery.</p>
<p><b>Power frequency (50/60 Hz) magnetic field</b> IEC 61000-4-8</p>	<p>30 A/m</p>	<p>30 A/m</p>	<p>Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.</p>
<p><b>Note:</b> UT is the a.c. mains voltage prior to application of the test level.</p>			

declaration - electromagnetic immunity

Cefar TENS is intended for use in the electromagnetic environment specified below. The customer or the user of Cefar TENS should assure that it is used in such an environment.

Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic environment - Guidance
<p><b>Conducted RF</b> IEC 61000-4-6</p>	<p>3 V 0.15 MHz to 80 MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz</p>	<p>3 V 0.15 MHz to 80 MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of Cefar TENS, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance</p> $d = \left[ \frac{3.5}{\sqrt{f}} \right] \sqrt{P}$ <p>150 KHz to 80 MHz</p> $d = \left[ \frac{3.5}{\sqrt{f}} \right] \sqrt{P}$ <p>80 MHz to 800 MHz</p> $d = \left[ \frac{7}{\sqrt{f}} \right] \sqrt{P}$ <p>80 MHz to 2.7 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

<b>Radiated RF</b> <b>IEC 61000-4-3</b>	10V/m 80 MHz to 2.7 GHz	10V/m	
<p><b>Note 1</b> At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p><b>Note 2</b> These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p><b>a</b> Field strengths from fixed RF transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Cefar TENS is used exceeds the applicable RF compliance level above, Cefar TENS should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating Cefar TENS.</p> <p><b>b</b> Over the frequency range 0.15 MHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

Recommended separation distances between portable and mobile RF communications equipment and Cefar TENS

Cefar TENS is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of Cefar TENS can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and Cefar TENS, as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	0.15 MHz to 80 MHz $d=1.2 \sqrt{P}$	80 MHz to 800 MHz $d=1.2 \sqrt{P}$	80 MHz to 2.7 GHz $d=2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer

**Note 1** At 80 MHz and 800 MHz, the higher frequency range applies.

**Note 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

declaration - IMMUNITY to proximity fields from RF wireless communications equipment

Cefar TENS is intended for use in an electromagnetic environment in which RF wireless communications equipment are controlled.

Immunity Test	IEC60601 test level				Compliance level	Electromagnetic environment - guidance
	Test frequency	Modulation	Maximum power	Immunity level		
Radiated RF IEC 61000-4-3	385 MHz	**Pulse Modulation: 18Hz	1.8W	27 V/m	27 V/m	
	450 MHz	*FM+ 5Hz deviation: 1kHz sine	2 W	28 V/m	28 V/m	
	710 MHz 745 MHz 780 MHz	**Pulse Modulation: 217Hz	0.2 W	9 V/m	9 V/m	
	810 MHz 870 MHz 930 MHz	**Pulse Modulation: 18Hz	2 W	28 V/m	28 V/m	
	1720 MHz 1845 MHz 1970 MHz	**Pulse Modulation: 217Hz	2 W	28 V/m	28 V/m	
	2450 MHz	**Pulse Modulation: 217Hz	2 W	28 V/m	28 V/m	
	5240 MHz 5500 MHz 5785 MHz	**Pulse Modulation: 217Hz	0.2 W	9 V/m	9 V/m	

**Note \*** - As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

**Note\*\*** - The carrier shall be modulated using a 50 % duty cycle square wave signal.

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DJO Global, Inc
U.S.A.