 chattanooga®

Intelect® **Transport 2**

Combo and Ultrasound

User Manual

EN

FOREWORD	4	INTELECT® TRANSPORT 2 SET COMPONENTS	14
INTENDED USER PROFILE	4	COMBO SET INCLUDES:	14
INTENDED ENVIRONMENT FOR USE	4	ULTRASOUND SET INCLUDES:	14
INTENDED USE	4	HEAD TO CART FIXATION	14
PRECAUTIONARY INSTRUCTIONS	4	INSTALLING THE BATTERY PACK	14
DESCRIPTION OF DEVICE MARKINGS	5	CHARGING THE BATTERY PACK	15
ELECTROTHERAPY INDICATIONS	6	USING THE BATTERY PACK	15
INDICATIONS	6	SYSTEM SPECIFICATIONS AND DIMENSIONS	16
CONTRAINDICATIONS	6	POWER	16
ADDITIONAL PRECAUTIONS	6	ELECTRO STIMULATION SPECIFICATIONS	16
ADVERSE EFFECTS	7	ULTRASOUND SPECIFICATIONS	16
ULTRASOUND INDICATIONS	7	OUTPUT POWER	16
INDICATIONS	7	ULTRASOUND SPATIAL PATTERN	17
CONTRAINDICATIONS	7	GENERAL SYSTEM OPERATING AND STORAGE TEMPERATURE	17
ADDITIONAL PRECAUTIONS	7	WAVEFORMS	18
ADVERSE EFFECTS	7	WAVEFORM GUIDANCE	22
CAUTION	8	ELECTROTHERAPY PATIENT PREPARATION AND ELECTRODE PLACEMENT	27
WARNING	9	ULTRASOUND PATIENT PREPARATION	28
DANGER	10	STARTING, STOPPING, AND INTERRUPTING ELECTROTHERAPY	28
PRODUCT DESCRIPTION	11	STARTING, STOPPING, AND INTERRUPTING ULTRASOUND THERAPY	30
COMPONENTS	11	STARTING, STOPPING, AND INTERRUPTING COMBINATION THERAPY	31
HEAD	11	USING CLINICAL INDICATIONS	32
CART (OPTIONAL)	11	CREATING A USER PROTOCOL FOR ELECTROTHERAPY	33
BATTERY MODULE (OPTIONAL)	11	CREATING A USER PROTOCOL FOR ULTRASOUND THERAPY	34
ULTRASOUND APPLICATORS	11		
OPERATOR INTERFACE	12		

CREATING A USER PROTOCOL FOR COMBINATION THERAPY	34
USING PROTOCOLS	35
SYSTEM UTILITIES	35
ERROR CODES	39
REPLACEMENT AND OPTIONAL ACCESSORIES	40
STANDARD ACCESSORIES	40
OPTIONAL ACCESSORIES	40
CLEANING THE INTELECT® TRANSPORT 2	41
YEARLY CALIBRATION	41
DEVICE DISPOSAL	41
WARRANTY REPAIR/OUT OF WARRANTY REPAIR	41
WARRANTY	42
ELECTROMAGNETIC COMPATIBILITY (EMC) TABLES	43

FOREWORD

This manual is intended for users of Intellect® Transport 2. It contains general information on operation, precautionary practices, and maintenance.

To maximize use, efficiency, and the life of the system, please read this manual thoroughly and become familiar with the controls, as well as the accessories before operating the system. In addition to the above information, this manual contains care and installation instructions for the optional Cart and the users of the Intellect® Transport 2.

Before administering any treatment to a patient, the users of this equipment should read, understand, and follow the information contained in this manual for each mode of treatment available, as well as the indications, contraindications, cautions, warnings, and dangers. Consult other resources for additional information regarding the application of electrotherapy and ultrasound.

INTENDED USER PROFILE

The intended user of this device is a licensed clinical professional. The user should be able to:

- Read and understand the operator's manual, warnings, cautions and dangers.
- Sense auditory and visual signals.
- Read and understand indications and contraindications of the device

INTENDED ENVIRONMENT FOR USE

The device is intended to be operated in the clinical setting including chiropractic clinics, physical therapist clinics or other rehabilitation settings.

 <p>ETL CLASSIFIED C US Intertek 5023416</p>	<p>The presence of this label indicates the machine was certified by ETL with the Statement: COMFORMS TO AAMI STD ES 60601-1, IEC STD 60601-1-6, 60601-2-5, 60601-2-10 CERTIFIED TO CSA STD C22.2 NO. 60601-1, No. 60601-1-6, No. 60601-2-5, No. 60601-2-10</p>
---	---

INTENDED USE

The Intellect® Transport 2 devices comprise of a range of multi-modality (TENS; NMES, Ultrasound) therapies intended to be used by healthcare professionals using TENS, NMES and Therapeutic Ultrasound for the treatment of various musculoskeletal and skeletal muscle deficit disorders.

The Intellect® Transport 2 product range offers the following models:

The Intellect® Transport 2 Ultrasound device delivering Therapeutic Ultrasound

The Intellect® Transport 2 Combo Device delivering both Electrotherapy (TENS and NMES) and Therapeutic Ultrasound either simultaneously or independently

PRECAUTIONARY INSTRUCTIONS

The precautionary instructions found in this section and throughout this manual are indicated by specific symbols. Intended users need to understand these symbols and their definitions before operating this equipment. The definition of these symbols are as follows:

Text with a "CAUTION" indicator explains possible safety infractions that have potential to cause minor or moderate injury or damage to the equipment.

Text with a "WARNING" indicator explains possible safety infractions that will potentially cause serious injury and equipment damage.

Text with a "DANGER" indicator will explain possible safety infractions that are imminently hazardous situations that would result in death or serious injury.

NOTE: Throughout this manual, "NOTE" indicators provide helpful information regarding the particular area of function being described.

DESCRIPTION OF DEVICE MARKINGS

The markings on the unit are assurance of its conformity to the highest applicable standards of medical equipment safety and electromagnetic compatibility and conform to ISO 7010 and ISO15-223-1 One or more of the following markings may appear on the device:

Refer to Instructional Manual Booklet		Relative Humidity Range	
Warning, Caution, or Danger		Atmospheric Pressure Range	
Electrical Type BF Equipment		Test agency	
Electrical Type B Equipment		Alternating current	
Ultrasound		IP21	IP21
Stim		WEEE Directive conformity	
Clinical library		Shelf life	
Increase for time or intensity	+	Batch number	LOT
Decrease for time or intensity	-	US amplitude modulated	
ON/OFF		MD	MD
Manufacturer		Start	
Date and country of manufacture		Stop	
Catalogue number	REF	Pause	
Serial number	SN	Parameter Display/Enter	
Fragile, handle with care		Back	
This end up		Up Arrow	
Keep dry		Down Arrow	
Temperature Range		Battery Indicator	
CE Mark of Conformity with notified body number	CE 2797	Charge Indicator	

ELECTROTHERAPY INDICATIONS

INDICATIONS

TENS:

- Symptomatic relief of acute and chronic pain associated with musculoskeletal conditions
- Management and alleviation of post-operative pain

NMES:

- Skeletal muscle deficit disorders resulting in benefits such as:
 - » Muscle re-education
 - » Maintaining/increasing range of motion

CONTRAINDICATIONS

The Intellect® Transport 2 should NOT be used under the following conditions:

- Do not use for symptomatic local pain relief unless etiology is established or unless a pain syndrome has been diagnosed.
- Do not use when cancerous lesions are present in the treatment area.
- Do not apply stimulation over swollen, infected, inflamed areas or skin eruptions (e.g., phlebitis, thrombophlebitis, varicose veins, etc.).
- Do not use when patient is suspected or known to have an infectious disease and/or disease where it is advisable, for general medical purposes, to suppress heat or fevers.
- Do not place electrode placements to the carotid sinus region (anterior neck) or transcranially (through the head).
- Do not use on pregnant women. Safety has not been established for the use of therapeutic electrical stimulation during pregnancy.
- Do not use Intellect® Transport 2 on patients who have or have had implantable neurostimulating cardiac demand pacemakers, ICD, or other implantable electronic devices.
- Do not use Intellect® Transport 2 on patients with body worn electromechanical medical devices, i.e. insulin pump.

- Do not use this system in an MRI or CT environment. The Intellect® Transport 2, its components, and accessories are not to be present in an MRI or CT environment.
- Do not apply stimulation transthoracically or on the chest, the introduction of electrical current into the heart may cause cardiac arrhythmia.
- Do not apply stimulation over the anterior neck or mouth. Severe spasm of the laryngeal and pharyngeal muscles may occur and the contractions may be strong enough to close the airway or cause difficulty in breathing.

ADDITIONAL PRECAUTIONS

- Use caution for patients with suspected or diagnosed heart problems.
- Use caution for patients with suspected or diagnosed epilepsy.
- Use caution in the presence of the following:
 - » When there is a tendency to hemorrhage following acute trauma or fracture
 - » Following recent surgical procedures when muscle contraction may disrupt the healing process
 - » Over a menstruating or pregnant uterus
 - » Over areas of the skin that lack normal sensation
- Powered muscle stimulators should be used only with the lead wires and electrodes recommended for use by the manufacturer.
- With TENS waveforms, isolated cases of skin irritation may occur at the site of electrode placement following long-term application.
- The effective management of pain by TENS waveforms is highly dependent upon patient selection by a person qualified in pain management.

ELECTROTHERAPY INDICATIONS (CONTINUED)

ADVERSE EFFECTS

- Skin irritation and burns beneath the electrodes have been reported with the use of powered muscle stimulators. The irritation can usually be reduced by using an alternative conductive medium or an alternative electrode placement
- Potential adverse effects with TENS are skin irritation and electrode burns

Note: Skin irritation and burns beneath the electrodes can be reduced or avoided by using appropriate electrode size and ensuring optimal contact quality.

Some people, with very sensitive skin, may experience redness under the electrodes after a session. Generally, this redness is totally harmless and usually disappears after 10 to 20 minutes. However, do not start another stimulation session on the same area if the redness is still visible.

- Do not use over the thoracic area if the patient is using a cardiac pacemaker.
- Do not use over a healing fracture.
- Do not use over or applied to the eye.
- Do not use over a pregnant uterus.
- Tissue necrosis might result if the device is used on ischemic tissues in individuals with vascular disease, where the blood supply would not keep up with the metabolic demand.
- Do not use Intelect® Transport 2 on patients who have or have had implantable neurostimulating cardiac demand pacemakers, ICD, or other implantable electronic devices.
- Do not use Intelect® Transport 2 on patients with body worn electro mechanical medical devices, i.e. insulin pump.
- Do not use this system in an MRI or CT environment. The Intelect® Transport 2, its components, and accessories are not to be present in an MRI or CT environment.

ULTRASOUND INDICATIONS

INDICATIONS

- Musculoskeletal conditions; facilitating pain relief to the affected area

CONTRAINDICATIONS

- Do not use for symptomatic local pain relief unless etiology is established or unless a pain syndrome has been diagnosed.
- Do not use when cancerous lesions are present in the treatment area.
- Do not use when patient is suspected or known to have infectious disease and/or disease where it is advisable, for general medical purposes, to suppress heat or fevers.
- Do not use over or near bone growth centers until bone growth is complete.

ADDITIONAL PRECAUTIONS

Additional precautions should be used when ultrasound is used on patients with the following conditions:

- Over an area of the spinal cord following a laminectomy, i.e., when major covering tissues have been removed
- Over anesthetized areas
- On patients with hemorrhagic diatheses

ADVERSE EFFECTS

Stinging sensation and temporary local sensitivity have been reported during or following ultrasound therapy.

 **CAUTION**

- Read, understand, and practice the precautionary and operating instructions. Know the limitations and hazards associated with using any electrical stimulation or ultrasound device. Observe the precautionary and operational decals placed on the unit
- The unit should be routinely checked before each use to determine that all controls function normally; especially that the intensity control properly adjusts the intensity of the electrotherapy and ultrasonic power output in a stable manner. Also, determine that the treatment time control actually terminates electrotherapy and ultrasonic power output when the timer reaches zero.
- This unit should be operated in temperatures between 41°F and 104°F (5°C and 40°C), with relative humidity ranging from 15%-90%, and where the atmospheric pressure is between 70 k Pa and 106 kPa.
- This unit should be transported and stored in temperatures between 41°F and 104°F (-20°C and 60°C), with relative humidity ranging from 10%-90%, and where the atmospheric pressure is between 50 kPa and 106 kPa.
- DO NOT disassemble, modify, or remodel the unit or accessories. This may cause unit damage, malfunction, electrical shock, fire, or personal injury.
- DO NOT use sharp objects such as a pencil point or ballpoint pen to operate the buttons on the control panel.
- DO NOT permit any foreign materials including, but not limited to, liquids such as water and cleaning agents, inflammables, and metallic objects to enter the unit. To prevent unit damage, malfunction, electrical shock, fire, or personal injury.
- DO NOT remove cover. There are no user-serviceable parts inside the unit. If a malfunction occurs, discontinue use immediately and consult the dealer for repair service.
- Before each use inspect cables and connectors.
- Handle the applicator with care. Inappropriate handling of the applicator may adversely affect its characteristics.
- Do not drop the applicator on hard surfaces. This is likely to damage the ultrasound head crystal. Damage resulting from this condition is not covered under the warranty.
- Before each use, inspect the applicator for cracks, which may allow the ingress of conductive fluid.
- Electrotherapy output current density is inversely related to electrode size. Always exercise caution with current densities more than 2mA/cm². Improper application may result in patient injury. If any question arises as to the proper electrode size, consult a licensed practitioner prior to therapy session.
- For waveforms with a DC component:
 - » Do not shave electrodes application area
 - » Warn the patient that tingling sensation under electrodes is normal and it is not linked to burn risk.
 - » Rinse thoroughly treatment area with tap water immediately after the treatment
- Where the integrity of the external protective earth conductor arrangement is in doubt, equipment shall be operated from its internal electrical power source, the battery.
- The battery pack should be removed when storing the unit for extended periods of time.
- Using a high intensity electrotherapy setting in conjunction with high intensity ultrasound setting may cause the unit to reset.
- Failure to use and maintain the Intellect Transport 2 device and its accessories in accordance with the instructions outlined in this manual will invalidate your warranty.

 **WARNING**

- Make certain the unit is electrically grounded by connecting only to a grounded electrical service receptacle conforming to the applicable national and local electrical codes.
- Transport 2 devices are designed to comply with electromagnetic safety standards. However, Transport 2 devices generate, use, and can radiate radio frequency energy and, if not installed and used in accordance with instructions for use, may cause harmful interference to other devices in the vicinity. Harmful interference to other devices can be determined by turning this equipment on and off. Try to correct the interference using one or more of the following:
 - » Reorient or relocate the receiving device
 - » Increase the separation between the equipment
 - » Connect the equipment to an outlet on a different circuit from that to which the other device(s) are connected
 - » Consult your authorized DJO dealer for help
- Do not use Transport 2 Combo device in conjunction with electronic monitoring equipment (such as ECG monitors and ECG alarms). Electronic monitoring equipment may not operate properly when electrical stimulation is in use.
- DO NOT operate Transport 2 in an environment where other devices are being used that intentionally radiate electromagnetic energy in an unshielded manner. For example
 - » DO NOT operate the Intelect® Transport 2 within the vicinity or environment of an ultrasonic diathermy system.
 - » DO NOT operate the Intelect® Transport 2 within the vicinity or environment of any microware and RF shortwave diathermy system.
- The energy from above systems can be transferred to patients with an implanted neurostimulation device causing tissue damage and resulting in severe injury or death even if the implanted neurostimulation system is turned “off.”
- Portable RF communications equipment should be used no closer than 30 cm (12 inches) to any part of the Intelect Transport 2, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- Simultaneous connection of a PATIENT to a high frequency surgical ME EQUIPMENT may result in burns at the site of the STIMULATOR electrodes and possible damage to the STIMULATOR.
- DO NOT operate the Intelect Transport 2 unit when connected to any unit other than DJO devices or accessories specifically described in this IFU as part of Transport 2 system or that have been specified as being compatible with the Transport 2. The use of other companies’ accessories, transducers or cables may result in increased emissions or decreased immunity of Transport 2 devices and resulting improper operations. DJO, LLC is not responsible for any consequence resulting from using products manufactured by other companies.
- DO NOT apply stimulation transthoracically because the introduction of electrical current into the heart may cause cardiac arrhythmia.
- DO NOT use contaminated, electrodes, lead wires, and gel which can lead to infection.
- DO NOT use electrode on multiple patients can lead to infection.
- Keep electrodes separated during treatment. Electrodes in contact with each other could result in improper stimulation or skin burns.
- DO NOT apply electro stimulation treatment during bath, shower, sauna.

WARNINGS (CONTINUED)

- The Intellect® Transport 2 may be susceptible to Electro-Static Discharge (ESD) at greater than ± 6 kV when first grasping the Ultrasound applicator. In the event of such a discharge, the Intellect® Transport 2 may display a permanent error. The Intellect® Transport 2 will terminate all active outputs (stim, ultrasound,), automatically place the unit in a safe state. Do not turn the unit on or off while it is connected to the patient.
- To prevent Electro-Static Discharge (ESD) at greater than ± 6 kV:
 - » Grasp and hold the Ultrasound prior to starting treatment. If the applicator must be put down prior to completion of treatment, stop the current treatment first and then place the applicator in the holder.
 - » Maintain humidity in the use environment to at least 50% relative humidity.
 - » Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, DJO recommends implementing additional controls to maintain relative humidity to at least 50%.
 - » Communicate these ESD-precautionary procedures to healthcare staff, contractors, visitors, and patients.
- In the event that an Error message or Warning appears beginning with a 2 or 3, immediately stop all use of the unit and contact the dealer or DJO, LLC for service. Errors and Warnings in these categories indicate an internal problem with the unit that must be tested by DJO, LLC or a Field Service Technician certified by DJO, LLC before any further operation or use of the system.
 - » Use of a unit that indicates an Error or Warning in these categories may pose a risk of injury to the patient, user, or extensive internal damage to the system.
- Disconnect the system from the power source and battery before attempting any maintenance, installation, removal, or replacement procedures to prevent electrical shock and possible damage to system.

Serious incidents occurring in relation to the device should be reported to the manufacturer and the local competent authority of the Member State in which the user and/or patient is established.

**DANGER**

- Device is not designed to be used in oxygen rich environment, Explosion hazard if the device is used in the presence of flammable anesthetic mixture with air, oxygen, or nitrous oxide.
- Do not reverse the polarity of the battery pack. Doing so can increase the individual cell temperature and cause cell rupture or leakage.

PRODUCT DESCRIPTION

The Intellect® Transport 2 comes in two different models: an ultrasound only unit and a combo unit with two channel electrotherapy and ultrasound.

This equipment is to be used only under the prescription and supervision of a licensed practitioner.

COMPONENTS

Throughout these instructions the terms "left" and "right" referring to the machine sides are from the perspective of a user standing in front of the unit.

The components of the Intellect Transport 2 (shown below) can be used on both the combo model and ultrasound model.

HEAD

CART (OPTIONAL)

BATTERY MODULE (OPTIONAL)

Battery is an 18V 3350mAh Li-Ion rechargeable battery.

ULTRASOUND APPLICATORS

Ultrasound applicators are available in 2 cm², 5 cm² and 10 cm² options.

OPERATOR INTERFACE

The Intellect Transport 2 Ultrasound device has two ultrasound connectors and software interface with only ultrasound capabilities. In addition to all of the ultrasound capabilities, combo devices have electrostimulation and electrostimulation-ultrasound combination features, as well as two connectors for electrostimulation lead wires.

User Interface

1. "ON/OFF" button. Controls the flow of electricity to the unit
Note: Make certain there are no electrodes on the patient when turning the unit on or off.
2. CLINICAL LIBRARY button
Select this button to access the following functions:
 - Retrieve User Protocol
 - Restore Factory Settings
 - Restore Factory Protocols
 - Language
 - View Unit Information
 - Backlight
3. TIME button
 - Press the Up or Down arrow buttons to set total treatment time of therapy
4. INTENSITY button
 - Use the up or down arrow on the INTENSITY button to increase or decrease output power.
5. 4-WAY button
 - Up Arrow
When the window displays a list of options, press the Up Arrow button to scroll up the list.
 - Down Arrow
When the window displays a list of options, press the Down Arrow button to scroll down the list.
 - Back
Use this button to return to the previous window.
 - Parameter Display/Enter
Select this button to display the parameters of the waveform during treatment. Also, this button is used to accept the highlighted selection.
6. STOP button
 - Select this button to stop a treatment session
7. START button
 - Select Start to begin a treatment session
8. PAUSE button
Use this button to pause the treatment session. To restart therapy, press the PAUSE button.
9. Ultrasound Applicator holder
10. Battery cover
11. Device handle

INTELECT® TRANSPORT 2 SET COMPONENTS

The components of the Intellect® Transport 2 set are shown below.

COMBO SET INCLUDES:

47544	INTELECT TRANSPORT 2 COMBO
47979	User Manual
47383	5 cm ² Ultrasound Applicator
48100	International Accessory Kit
- 40-0097	Leadwire Assembly, Ch 1
- 40-0098	Leadwire Assembly, Ch 2
- 79967-4	Silicon/Carbon Electrode 6x8cm Black
- 79970-4	Sponge Pocket 9x10.5 cm
- 12-10648	Nylatex Wrap 2-1/2x24 Sewn
- 42197	Dura-Stick+ 5cm Round Electrodes

ULTRASOUND SET INCLUDES:

47242	Intellect Transport 2 Ultrasound
47979	User Manual
47383	5 cm ² Ultrasound Applicator

HEAD TO CART FIXATION

The optional Therapy System Cart, is designed for use with the Intellect® Transport 2 device and allows the user to easily transport the system from patient to patient within the clinic as well as store all necessary accessories, supplies, and applicators used for the various modalities of the System. The fixation of the head to the cart is magnetic.

Remove the Intellect® Transport 2 device and cart from the shipping carton. Visually inspect for damage. Report any damage to the carrier immediately.

To assemble the Transport 2 Head to the Cart, follow these steps:

1. Insert the device onto the top of the cart by sliding the front bottom of the device into the plastic lip on the top front of the cart.
2. Release device back gently on the cart. Magnets will help to position the device correctly on the cart top.

INSTALLING THE BATTERY PACK

The Intellect Transport 2 devices accommodate both AC mains power and an optional Li-ion battery pack. To install the battery pack in a Intellect Transport 2 device, do the following:

1. Locate the battery access door at the bottom of the unit and loosen the screws with a coin or a flat head screw driver.

2. Remove the battery access door and retain it.
3. Connect the battery pack cable to the unit's battery connector in the bottom of the battery recess. Put the battery pack into the unit.
4. Replace the battery access door and re-tighten the screws using the coin or the screwdriver.
5. Reverse the steps in this section in order to remove the battery pack.

CHARGING THE BATTERY PACK

The battery pack is automatically charged by the unit whenever there is mains power connected. Charging may be interrupted during operation of the unit by the control circuitry to limit total power consumption. A fully charged battery will provide 1.5-3 hours of treatment depending on the applicator and the pulsed mode used.

NOTE: Even when the battery pack is connected, the unit will default to mains power when plugged in.

USING THE BATTERY PACK

To save battery power, the Intelect Transport 2 devices are equipped with a "power off" function. This function is activated when the unit is powered on and has been left idle for approximately 5 minutes, at which time the unit powers off. To restore power, press the Power On/Off button.

SYSTEM SPECIFICATIONS AND DIMENSIONS

	Width	Depth	Height	Weight (no battery)
Intelect Transport 2 Head Unit				
COMBO	25.0 cm	34.6 cm	11.9 cm	2.2 kg
UltraSound	25.0 cm	34.6 cm	11.9 cm	1.9 kg
Cart configurations				
Cart (Safe working load 6.5 kg)	48 cm (MAX)	52 cm (MAX)	96 cm	10.1 kg
Device on cart	-	-	108.5 cm	-

POWER

Input	100 - 240 V AC, 1.0 to 0.42 A, 50/60 Hz 100 W Max
Electrical Class	CLASS I
Mode of Operation	Continuous

Electrical Type (Degree of Protection)

Ultrasound	.TYPE B
Electrotherapy	.TYPE BF

ELECTRO STIMULATION SPECIFICATIONS

Output specifications are described for each waveform from pages 23-25.

Unless otherwise specified, electrotherapy controls accuracy is: $\pm 20\%$.

Load impedance: 500-1000 Ohm

CC = constant current, effect of load impedance on voltage

CV = constant voltage, effect of load impedance on current

ULTRASOUND SPECIFICATIONS

Frequency	1 MHz; 3.3 MHz
Duty Cycles	10%, 20%, 50%, Continuous
Pulse Repetition Rate	16, 48, 100 Hz
Pulse duration:	1 -5 ms
	Max (ON): 5 ms
	Min (OFF): 5 ms

OUTPUT POWER


US applicator Frequency	2 cm ²		5 cm ²		10 cm ²	
	1 MHz	3.3 MHz	1 MHz	3.3 MHz	1 MHz	3.3 MHz
Effective Radiating Area ERA INTL (cm ²)	1.5	1.5	3.0	3.0	5.8	5.8
Max Output power in Continuous mode	3W	3W	6W	6W	12.7W	6.9W
Max Output power in Pulsed mode	4.5W	4.5W	9W	9W	17.4W	13.9W
Max Intensity in Continuous mode	2W/cm ²	2W/cm ²	2W/cm ²	2W/cm ²	2.2W/cm ²	1.2W/cm ²
Max Intensity in Pulsed mode	3W/cm ²	3W/cm ²	3W/cm ²	3W/cm ²	3W/cm ²	2.4W/cm ²

*) Max output power has tolerance of $\pm 20\%$ OR $\pm 0.25[W]$ whichever is greater

Unless otherwise specified, ultrasound controls accuracy is:	$\pm 20\%$
Peak to Average Ratio:	1:1, at 50% Duty Cycle 4:1, at 20% Duty Cycle 9:1, at 10% Duty Cycle
Beam Nonuniformity Ratio	<6:1
Beam Type	Collimating
Treatment Time	1 to 30 min

ULTRASOUND SPATIAL PATTERN

The following chart represents the distribution of the ultrasonic radiation field and the orientation of the field with respect to the 5cm² applicator (Y-plane represents voltage in Vrms and X-plane represents applicator head surface in 1mm resolution).

IPXX Rating for Unit

Rated to IP21

IP2* Protection against fingers or other object not greater than 80mm in length and 12mm in diameter

*1 Protection from vertically dripping water

IPXX Rating for US applicator

Rated to IPX7

IPX7 Protection from immersed in water (up to 1m depth)

GENERAL SYSTEM OPERATING AND STORAGE TEMPERATURE

Operating Conditions

The device will meet its requirement under the following conditions:

Temperature:	5°C to 40°C
Relative Humidity:	15% to 90%
Atmospheric Pressure:	70 kPa to 106 kPa

Transport and Storage Conditions

The device will remain in proper condition under the following conditions:

Temperature:	-20°C to 60°C
Relative Humidity:	0% to 90%
Atmospheric Pressure:	50 kPa to 106 kPa

Time required for the Intellect[®] Transport 2 to warm from the minimum storage temperature between uses until the Intellect[®] Transport 2 is ready for its INTENDED USE when the ambient temperature is 20°C: 5H

WAVEFORMS

CC: Constant Current

CV: Constant Voltage

The software includes a Dynamic Output feature to ensure that the electrotherapy output to the patient adheres to the limits specified in IEC 60601-2-10:2012 section 201.12.4.104 which specifies the maximum RMS current allowed to the patient.

The software achieves this limitation by performing the following whenever the user changes the desired output intensity:

. The software will calculate the RMS current that would be delivered to the patient using the newly set parameters (i.e. desired output intensity, phase duration, frequency, etc.).

. If this calculated RMS current value exceeds the maximum allowed for the current waveform (defined in IEC 60601-2-10:2012 section 201.12.4.104 Table 201.101) then the user will not be allowed to set the intensity to this new value.

WARNING

Never use electrodes smaller than 2" when applying IFC treatment to a patient.

IFC (INTERFERENTIAL) TRADITIONAL (4 POLE)

Interferential Current is a medium frequency waveform. Current is distributed through two channels (four electrodes). The currents cross each other in the body at the area requiring treatment. The two currents interfere with each other at this crossing point, resulting in a modulation of the intensity (the current intensity increases and decreases at a regular frequency).

Output Mode	Electrodes
Available on Channel	1, 2
Treatment Time	1-60 min
Mode Selection	CC
Output Intensity	0-100 mA (CC) 0-100 V (CV)
Beat Frequency	1-200 Hz
Carrier Frequency	2000-10,000 Hz
Cycle Time	Continuous or User Defined
Sweep Time	14 sec
Sweep Low Beat Frequency	1-199 Hz
Sweep High Beat Frequency	2-200 Hz
Scan Percentage	Static, 40%, 100%, Manual
IRMS	0-78mA
DC component	No

IFC PREMODULATED (TRADITIONAL 2 POLE)

Premodulated Current is a medium frequency waveform. Current comes out of one channel (two electrodes). The current intensity is modulated: it increases and decreases at a regular frequency (the Amplitude Modulation Frequency).

Output Mode	Electrodes
Available on Channel	1, 2
Treatment Time (STIM)	1-60 min
Treatment Time (COMBO)	1-30 min
Mode Selection	CC or CV
Output Intensity	0-100 mA (CC) 0-100 V (CV) Carrier
Beat Fixed (Sweep Off)	1-200 Hz
Cycle Time	Continuous or User Defined
Frequency	2,000-10,000 Hz
Sweep Low Beat Frequency	1- 199 Hz
Sweep High Beat Frequency	2-200 Hz
IRMS	0-55mA
DC component	No

RUSSIAN

Russian Current is a sinusoidal waveform, delivered in bursts or series of pulses.

Output Mode	Electrodes
Available on Channels	1, 2
Treatment Time	1-60 min
Mode Selection	CC or CV
Output Intensity	0-100 mA (CC) 0-100 V (CV)
Burst Frequency	1-100 bps
Carrier Frequency	2,500 Hz
Cycle Time	Continuous or User Defined
Duty Cycle	10%, 20%, 30%, 40%, 50%
Ramp	0-5 sec
IRMS	0-39mA
DC component	No

WAVEFORMS (CONTINUED)

VMS™ is a symmetrical biphasic waveform with a 100 sec interphase interval. Because the pulse is relatively short, the waveform has a low skin load, making it suitable for applications requiring high intensities, such as in muscle re-education protocols.

Output Mode	Electrodes
Available on Channels	1, 2
Treatment Time (Stim)	1-60 min
Treatment time (Combo)	1-30 min
Mode Selection	CC or CV
Output Intensity	0- 140 mA (CC) 0-140 V (CV)
Anti-Fatigue	Off or On
Channel Mode	Single, Reciprocal, Co-Contract
Cycle Time	Continuous or User Defined
Frequency	1-200 pps
Phase Duration	30-1,000 sec
Ramp	0-5 sec
Set Intensity	Individual/both Channel Intensity Setting in Reciprocal and Co-Contract modes

Microcurrent is a monophasic waveform of very low intensity.

Output Mode	Electrodes
Available on channels	1, 2
Treatment Time	1-60 min
Mode Selection	CC
Output Intensity	0-1,000 A
Duty Cycle	50%
Frequency	0.1-1,000 pps
Polarity	Positive, Negative, or Alternating
IRMS	0-1mA
DC component	No

HIGH VOLTAGE PULSED CURRENT (HVPC)

The High Voltage Pulsed Current (HVPC) has a very brief pulse duration characterized by two distinct peaks delivered at high voltage. The waveform is monophasic (current flows in one direction only). The high voltage causes a decreased skin resistance making the current comfortable and easy to tolerate.

Output Mode	Electrodes
Available on Channels	1, 2
Treatment Time (Stim)	1-60 min
Treatment Time (Combo)	1-30 min
Mode Selection	CV
Output Intensity	0-500 V (CV)
Cycle Time	Continuous or User Defined
Display	Volts
Frequency	1-200 pps
Polarity	Positive or Negative
Ramp	0.5-5 sec
Sweep time	14sec
Sweep High Frequency	2-200 pps
Sweep Low Frequency	1-199 pps
IRMS	0-45mA
DC component	0-1.5mA

TENS- ASYMMETRICAL BIPHASIC

The Asymmetrical Biphasic waveform has a short pulse duration. It is capable of strong stimulation of the nerve fibers in the skin as well as of muscle tissue. This waveform is often used in TENS devices. Because of its short pulse, the patient typically tolerates the current well, even at relatively high intensities.

Output Mode	Electrodes
Output Intensity	0-140 mA (CC) 0-140 V (CV)
Available on Channel	1, 2
Treatment Time (Stim)	1-60 min
Treatment Time (Combo)	1-30 min
Mode Selection (Stim)	CC or CV
Mode Selection (Combo)	CV
Amplitude Modulation	0% (off) to 100% on 10% steps
Burst Frequency	0-10 bps
Cycle Time	Continuous or User Defined
Frequency	1-200 pps
FrequencySweep	On/Off
Phase Duration	30-400 sec
Sweep time	14 sec
Sweep Low Frequency	1-199 pps
Sweep High Frequency	2-200 pps
IRMS	0-50mA
DC component	No

WAVEFORMS (CONTINUED)

TENS- SYMMETRICAL BIPHASIC

The Symmetrical Biphasic waveform has a short pulse duration and is capable of strong stimulation of nerve fibers in the skin and in muscle. This waveform is often used in portable muscle stimulation units, and some TENS devices.

Output Mode	Electrodes
Available on Channel	1, 2
Treatment Time (Stim)	1-60 min
Treatment Time (Combo)	1-30 min
Mode Selection (Stim)	CC or CV
Mode Selection (Combo)	CV
Output Intensity	0-140 mA (CC) 0-140 V (CV)
Amplitude Modulation	0% (off) to 100% on 10% steps
Burst Frequency	0-10 bps
Cycle Time	Continuous or User Defined
Frequency	1-200 pps
Frequency Sweep	On/Off
Phase Duration	30-400 sec
Ramp	0-5 sec
Sweep Time	14sec
Sweep Low Frequency	1-199 pps
Sweep High Frequency	2-200 pps
IRMS	0-50 mA
DC component	No

GALVANIC: CONTINUOUS

Galvanic Current is a direct current flowing in one direction only. The current can be continuous or interrupted.

Output Mode	Electrodes
Available on Channels	1, 2
Treatment Time	1-60 min
Mode Selection	CC
Output Intensity	0-40 mA (CC)
Cycle Time	Continuous, or User Defined
Polarity Reversal	On or Off
	With Polarity Reversal On, Polarity will change in the middle of the treatment time.
IRMS	0-44 mA
DC component	Yes

TRÄBERT (ULTRAREIZ)

The Träbert Current is a monophasic waveform with a phase duration of 2 ms and a pause of 5 ms resulting in a frequency of approximately 143 Hz.

Output Mode	Electrodes
Available on Channels	1, 2
Treatment Time	1-60 min
Mode Selection	CC
Output Intensity	0-80 mA (CC)
Frequency	143 pps
Polarity Reversal	On or Off
	With Polarity Reversal On, Polarity will change in the middle of the treatment time.
Phase Duration	2 ms
IRMS	0-47mA
DC component	Yes

GALVANIC: INTERRUPTED

Galvanic Current is a direct current flowing in one direction only. The current can be continuous or interrupted.

Output Mode	Electrodes
Available on Channels	1, 2
Treatment Time	1-60 min
Mode Selection	CC
Output Intensity	0-40 mA (CC)
Pulse Duration	136 sec
Phase Interval	25 usec
Polarity Reversal	On or Off
	With Polarity Reversal On, Polarity will change in the middle of the treatment time.
Polarity Reversal Ramp	1 sec
IRMS	0-41 mA
DC component	Yes

WAVEFORMS (CONTINUED)

MONOPHASIC:

MONOPHASIC TRIANGULAR PULSED

The Monophasic Triangular Pulsed waveform is an interrupted unidirectional current with a triangular pulse shape.

Output Mode	Electrodes
Available on Channels	1, 2
Treatment Time	1-60 min
Mode Selection	CC
Output Intensity	0-60 mA (CC)
Phase Duration	0.1-1,000 ms
Phase Interval	5-5,000 ms
IRMS	0-27 mA
DC component	Yes

MONOPHASIC:

MONOPHASIC RECTANGULAR PULSED

The Monophasic Rectangular Pulsed waveform is an interrupted unidirectional current with a rectangular pulse shape.

Output Mode	Electrodes
Available on Channels	1, 2
Treatment Time	1-60 min
Mode Selection	CC
Output Intensity	0-60 mA (CC)
Phase Duration	0.1-1,000 ms
Phase Interval	5-5,000 ms
IRMS	0-47 mA
DC component	Yes

DIADYNAMIC WAVEFORMS

The Diadynamic waveforms are rectified alternating currents. The alternating current is modified (rectified) to allow the current to flow in one direction only.

Output Mode	Electrodes
Available on channels	1, 2
Treatment Time	1-60 min
Mode Selection	CC
Output Intensity	0-60 mA
MF (Monophasé Fixe) - Frequency of 50 Hz: phase duration of 10 ms followed by a pause of 10 ms.	
IRMS [mA]	0-33 mA
DF - Frequency of 100 Hz: phase duration of 10 ms	
CP - 1 second of MF followed abruptly by 1 second of DF.	
LP - Rhythmical fluctuation between 2 MF currents.	
CP-iso - A combination of MF and DF waveforms. CP-id: Same as CP-iso.	
IRMS	0-47 mA
DC component	Yes

WAVEFORM GUIDANCE

Table 1: Top-5 Recommended waveforms for each indication, with a ranking for the grade of recommendation

Ranking	Symptomatic relief of acute pain associated with musculoskeletal conditions	Symptomatic relief of chronic pain associated with musculoskeletal conditions	Management and alleviation of post operative pain	Muscle re-education	Maintaining/increasing range of motion
1	TENS Asym BP	TENS Asym BP	TENS Asym B	VMS	VMS
2	TENS Sym BP	TENS Sym BP	TENS Sym BP	–	–
3	IFC-4p, IFC-2p	–	IFC-4p, IFC-2p	Russian current	Russian current
4	VMS	IFC-4p, IFC-2p	VMS	TENS Sym BP, TENS Asym BP	TENS Sym BP, TENS Asym BP
5	HVPC	VMS	HVPC	Pulsed MP Rect / Triang	Pulsed MP Rect / Triang

Table 2: Recommended parameters for each waveform in each of the indications

NO.	Waveform	Default settings	Symptomatic relief of acute pain associated with musculoskeletal conditions			
		Note: Default settings are the same for any conditions. Please adjust according to this table	Freq	Intensity	Duty cycle	Treatment Time
1	TENS-Asymmetrical Biphasic	80Hz; cont; 20 min	80-150Hz	strong sensation/ no pain	continuous	when in pain
2	TENS- Symmetrical Biphasic	80Hz; cont; 20 min	80-150Hz	strong sensation/ no pain	continuous	when in pain
3	VMS™	50Hz; 2s ramp; cont; 20 min	80-150Hz	strong sensation/ no pain	continuous	when in pain
4	IFC (Interferential) Traditional (4 Pole)	4000Hz; 80/150Hz; manual; cont; 20 min	use default settings and increase intensity to strong sensory level/no pain			when in pain
5	IFC Premodulated (Traditional 2 Pole)	2500Hz; 80/150Hz; cont; 20 min	use default settings and increase intensity to strong sensory level/no pain			when in pain
6	Russian	2500Hz; 50Hz; DC 50%; cont; 20 min				
7	High Voltage Pulsed Current (HVPC)	100Hz; cont; 20 min	use default settings and increase intensity to strong sensory level/no pain			
8	MONOPHASIC: Monophasic Rectangular pulsed	1.0ms/500ms; 5 min				
9	MONOPHASIC: Monophasic Triangular pulsed	1.0ms/500ms; 5 min				

10	Diadynamic Waveforms	Diadyn settings (fixed), 10 min	Use default Diadynamic settings and increase intensity to strong sensory level/no pain			
11	Träbert (Ultrareiz)	Träbert settings (fixed), 10 min	Use default Träbert settings and increase intensity strong to sensory level/no pain			
12	Microcurrent	Neg; 1Hz; 20 min	Use default settings and increase intensity to 40-100µA			when in pain/ multiple hours
13	GALVANIC: Interrupted	Pos; cont; 10 min				
14	GALVANIC: Interrupted	Pos; cont; 10 min				

NO.	Waveform	Default settings	Symptomatic relief of acute pain associated with musculoskeletal conditions			
		Note: Default settings are the same for any conditions. Please adjust according to this table	Freq	Intensity	Duty cycle	Treatment Time
1	TENS-Asymmetrical Biphasic	80Hz; cont; 20 min	2-5 Hz	muscle twitches	continuous	30 min
2	TENS- Symmetrical Biphasic	80Hz; cont; 20 min	2-5 Hz	muscle twitches	continuous	30 min
3	VMS™	50Hz; 2s ramp; cont; 20 min	2-5 Hz	muscle twitches	continuous	30 min
4	IFC (Interferential) Traditional (4 Pole)	4000Hz; 80/150Hz; manual; cont; 20 min	4000Hz; 2-5Hz	muscle twitches	continuous	30 min
5	IFC Premodulated (Traditional 2 Pole)	2500Hz; 80/150Hz; cont; 20 min	4000Hz; 2-5Hz	muscle twitches	continuous	30 min
6	Russian	2500Hz; 50Hz; DC 50%; cont; 20 min				
7	High Voltage Pulsed Current (HVPC)	100Hz; cont; 20 min	use default settings and increase intensity to strong sensory level/no pain			
8	MONOPHASIC: Monophasic Rectangular pulsed	1.0ms/500ms; 5 min				
9	MONOPHASIC: Monophasic Rectangular pulsed	1.0ms/500ms; 5 min				
10	Diadynamic Waveforms	Diadyn settings (fixed), 10 min	Use default Diadynamic settings and increase intensity to strong sensory level/no pain			
11	Träbert (Ultrareiz)	Träbert settings (fixed), 10 min				
12	Microcurrent	Neg; 1Hz; 20 min	Use default settings and increase intensity to 40-100µA			when in pain/ multiple hours

13	GALVANIC: Continuous	Pos; cont; 10 min	Use default settings, with polarity reversal and increase intensity to mild sensory level/no pain
14	GALVANIC: Interrupted	Pos; cont; 10 min	Use default settings and increase intensity to mild sensory level/no pain

NO.	Waveform	Default settings	Management and alleviation of post operative pain			
		Note: Default settings are the same for any conditions. Please adjust according to this table	Freq	Intensity	Duty cycle	Treatment Time
1	TENS- Asymmetrical Biphasic	80Hz; cont; 20 min	80-150Hz	strong sensation/ no pain	continuous	when in pain
2	TENS- Symmetrical Biphasic	80Hz; cont; 20 min	80-150Hz	strong sensation/ no pain	continuous	when in pain
3	VMS™	50Hz; 2s ramp; cont; 20 min	80-150Hz	strong sensation/ no pain	continuous	when in pain
4	IFC (Interferential) Traditional (4 Pole)	4000Hz; 80/150Hz; manual; cont; 20 min	Use default settings and increase intensity to strong sensory level/no pain			when in pain
5	IFC Premodulated (Traditional 2 Pole)	2500Hz; 80/150Hz; cont; 20 min	Use default settings and increase intensity to strong sensory level/no pain			when in pain
6	Russian	2500Hz; 50Hz; DC 50%; cont; 20 min				when in pain
7	High Voltage Pulsed Current (HVPC)	100Hz; cont; 20 min	Use default settings and increase intensity to strong sensory level/no pain			
8	MONOPHASIC: Monophasic Rectangular pulsed	1.0ms/500ms; 5 min				
9	MONOPHASIC: Monophasic Rectangular pulsed	1.0ms/500ms; 5 min				
10	Diadynamic Waveforms	Diadyn settings (fixed), 10 min	Use default Diadynamic settings and increase intensity to strong sensory level/no pain			
11	Träbert (Ultrareiz)	Träbert settings (fixed), 10 min	Use default Träbert settings and increase intensity strong to sensory level/no pain			
12	Microcurrent	Neg; 1Hz; 20 min	Use default settings and increase intensity to 40-100µA			when in pain/ multiple hours
13	GALVANIC: Continuous	Pos; cont; 10 min				
14	GALVANIC: Interrupted	Pos; cont; 10 min				

NO.	Waveform	Default settings	Muscle re-education			
		Note: Default settings are the same for any conditions. Please adjust according to this table	Freq	Intensity	Duty cycle	Treatment Time

1	TENS- Asymmetrical Biphasic	80Hz; cont; 20 min	35-80Hz	muscle contraction; max tolerable	1:5 to 1:1 with on-time 5s to 20s	10 - 15 min
2	TENS- Symmetrical Biphasic	80Hz; cont; 20 min	35-80Hz	muscle contraction; max tolerable	1:5 to 1:1 with on-time 5s to 20s	10 - 15 min
3	VMS™	50Hz; 2s ramp; cont; 20 min	35-80Hz	muscle contraction; max tolerable	1:5 to 1:1 with on-time 5s to 20s	10 - 15 min
4	IFC (Interferential) Traditional (4 Pole)	4000Hz; 80/150Hz; manual; cont; 20 min				
5	IFC Premodulated (Traditional 2 Pole)	2500Hz; 80/150Hz; cont; 20 min				
6	Russian	2500Hz; 50Hz; DC 50%; cont; 20 min	Use default settings and increase intensity to max tolerable motor level			
7	High Voltage Pulsed Current (HVPC)	100Hz; cont; 20 min				
8	MONOPHASIC: Monophasic Rectangular pulsed	1.0ms/500ms; 5 min	Use default settings and increase phase duration (denervated muscles) & intensity to obtain muscle contraction			
9	MONOPHASIC: Monophasic Rectangular pulsed	1.0ms/500ms; 5 min	Use default settings and increase phase duration (denervated muscles) & intensity to obtain muscle contraction			
10	Diadynamic Waveforms	Diadyn settings (fixed), 10 min				
11	Träbert (Ultrareiz)	Träbert settings (fixed), 10 min				
12	Microcurrent	Neg; 1Hz; 20 min				
13	GALVANIC: Continuous	Pos; cont; 10 min				
14	GALVANIC: Interrupted	Pos; cont; 10 min				


NO.	Waveform	Default settings	Maintaining/increasing range of motion			
		Note: Default settings are the same for any conditions. Please adjust according to this table	Freq	Intensity	Duty cycle	Treatment Time
1	TENS- Asymmetrical Biphasic	80Hz; cont; 20 min	35-50Hz	muscle contraction producing ROM	5s/5s	20 min
2	TENS- Symmetrical Biphasic	80Hz; cont; 20 min	35-50Hz	muscle contraction producing ROM	5s/5s	20 min
3	VMS™	4000Hz; 80/150Hz; manual; cont; 20 min	Use default settings and increase intensity to produce muscle contraction moving the joint through its ROM			
4	IFC (Interferential) Traditional (4 Pole)	4000Hz; 80/150Hz; manual; cont; 20 min				
5	IFC Premodulated (Traditional 2 Pole)	2500Hz; 80/150Hz; cont; 20 min				
6	Russian	2500Hz; 50Hz; DC 50%; cont; 20 min	Use default settings and increase phase duration (denervated muscles) & intensity to obtain muscle contraction			
7	High Voltage Pulsed Current (HVPC)	100Hz; cont; 20 min				
8	MONOPHASIC: Monophasic Rectangular pulsed	1.0ms/500ms; 5 min				
9	MONOPHASIC: Monophasic Rectangular pulsed	1.0ms/500ms; 5 min				
10	Diadynamic Waveforms	Diadyn settings (fixed), 10 min				
11	Träbert (Ultrareiz)	Träbert settings (fixed), 10 min				
12	Microcurrent	Neg; 1Hz; 20 min				
13	GALVANIC: Continuous	Pos; cont; 10 min				
14	GALVANIC: Interrupted	Pos; cont; 10 min				

ELECTROTHERAPY PATIENT PREPARATION AND ELECTRODE PLACEMENT

- Examine the skin for any wounds and clean the skin.
- Apply the electrodes to the treatment area.
- Ensure the electrodes are applied securely to the skin.
- Ensure good contact between each electrode and the skin.
- Check the electrode contact regularly during the treatment.
- Examine the skin again after the treatment.
- Choose electrodes that fit the anatomy.
- View the Electrode Placement recommendations in the Treatment Review screen for the particular modality being used for treatment as a reference point only prior to administering treatment.
- Follow electrode manufacturer instructions.
- Please note the smaller the electrode size the higher the current density.

DURA-STICK® Electrodes

DURA-STICK® Electrodes are a self adhesive, disposable product designed specifically for use with Intellect® Transport 2. It is recommended that DURA-STICK® Electrodes be used whenever possible to ensure the highest level of contact with the treatment area and most uniform delivery of the prescribed electrotherapy treatment.


For Electrotherapy operation refer to page 28.

DURA-STICK® Electrode Instructions

Connecting Lead Wires

1. Insert the lead with the Red (+) electrode connector into one DURA-STICK® Electrode.
2. Insert the lead with the Black (-) electrode connector into the other electrode.
3. Make certain the lead wires are seated completely into the electrodes.

NOTE: Use of conductive medium or sponges is not required or recommended. DURA-STICK® electrodes are manufactured to ensure the optimum conductivity during therapy when properly applied.

Securing Electrodes

1. Remove the DURA-STICK® Electrodes from the protective backing.
2. Apply to the treatment area as prescribed.
3. Ensure the entire electrode surface is in contact with patient skin by pressing into place.

ULTRASOUND PATIENT PREPARATION

1. Examine the skin for any wounds and clean the skin.

Applicator Preparation and Use

2. Clean applicator before each therapy session with warm soapy water, check the applicator has no cracks prior to use.
3. Liberally apply transmission gel to the treatment area on the patient.
4. Move the applicator during therapy session in a circular motion. The area treated should be:
 - Twice the diameter of the applicator
 - For 5cm² US applicator: three times the diameter of the applicator if output power > 4 W, Continuous mode.

The applicator should always be held by the grip and not by the Ultrasound Applicator head.

5. If US Coupling is "On", the Applicator is properly coupled to the patient and administering ultrasound when the LED is constantly illuminated. If the applicator head becomes uncoupled the LED on the head will flash. If "US coupling" setting is ON, several beeps will be also heard until the head is coupled again. Treatment time stops during uncoupling.

NOTE: Ultrasound output will continue to be emitted in all US coupling modes even if the applicator is uncoupled. The output power is reduced to a very low level to prevent ultrasound head warming.

STARTING, STOPPING, AND INTERRUPTING ELECTROTHERAPY

The operator is able to navigate between treatment options, view parameter options on the display and make selections by pressing the buttons on the control panel. The display will provide continuous information during the treatments concerning power and elapsed time. Parameters are adjusted using control panel buttons on the front of the unit. The output can be stopped by pressing the "PAUSE" or "STOP" buttons located on the control panel.

WARNING

Do not turn the unit on or off while it is connected to the patient.

Note: Electrotherapy related manual items are for Combo device only

To apply electrotherapy, do the following:

1. Follow all appropriate procedures listed in the section entitled "**Electrotherapy Patient Preparation**" beginning on page 27.
2. Turn system power "ON" by pressing the Power On/Off button. The message "Initializing System" displays while the unit goes through self diagnostics until the home screen appears.
3. Connect the Lead Wires to the appropriate electrodes. To see a list of recommended electrodes and their preparation, see **page 26**.


NOTE: Do not use unnecessary force to connect the electrodes to the lead wires.

4. Place the self adhesive electrodes on the sites prescribed by a qualified practitioner. Make sure you press them firmly on the patient's skin to ensure good conductivity.
5. Depending on the type of waveform you intend to use and the number of patients you intend to treat, insert the Lead Wire into Channel 1, Channel 2, or both Lead Wire Connections on the Accessory Panel.


STARTING, STOPPING, AND INTERRUPTING ELECTROTHERAPY (CONTINUED)

6. Use the Up and Down Arrow buttons to highlight either Stim Channel 1 or Stim Channel 2.
7. Press the Enter button.
The Waveform screen displays.
8. Use the Up and Down Arrow buttons to highlight the appropriate waveform.

9. Press the Enter button.
The Electrotherapy parameter screen displays.

10. Press the Up and Down Arrow buttons to highlight the parameter you want to change and adjust it accordingly, making sure to press the Enter button after each adjustment.
11. Press Enter to accept the parameters.
The Parameter screen refreshes and the new parameters are displayed.
12. To begin therapy, press the START button.
The timer counts down, the output power ramps up below the timer, and "Running" displays below next to the timer.
NOTE: When the therapy time has expired, the unit beeps three times.
During therapy, you can press the TIME button to raise or lower treatment time (in one minute increments) using the up and down arrows.
During therapy, you can press the INTENSITY button to raise or lower the output using the up and down arrows.
Therapy can be interrupted at any time by pressing the STOP or PAUSE buttons.

When the STOP button is pressed, the unit stops emitting output, and the unit returns to the home screen. To resume therapy, repeat **steps 6-11**.

During treatment, the following occurs whenever the PAUSE button is pressed:

- the timer pauses
- the unit beeps quickly 5 times
- "Paused" displays below the timer
- the unit stops emitting output

To resume therapy, press the PAUSE button or the START button.

13. When you have completed treatment, remove all electrodes from the patient.
14. Turn system power "OFF" by pressing the Power On/Off button. The unit beeps once and the blue light on the Power On/Off button flashes intermittently.

STARTING, STOPPING, AND INTERRUPTING ULTRASOUND THERAPY

To apply ultrasound therapy, do the following:

1. Ensure the ultrasound applicator is plugged in.
2. On the Home screen, use the Up and Down Arrow buttons to highlight Ultrasound.
3. Press the Enter button.


The Ultrasound parameters screen displays.

4. Press the Up and Down Arrow buttons to highlight Frequency.

5. Press the Enter button to select 1 or 3 MHz. When the button is pressed, the frequency will toggle from 1 to 3.3 MHz and back again as long as the button is being pressed (excluding the 1 cm² sound head).

An audible tone will be heard when changes are made.

NOTE: With 2 cm² and 5 cm² sound heads, switching from 1 to 3.3 MHz results in no change in power. When using a 10 cm² head with greater than 10 watts, changing from 1 to 3.3 MHz reduces power to 10 watts.

6. Press the Up and Down Arrow buttons to highlight Duty Cycle.
7. Press the Enter button to display the Duty Cycle menu.
8. Press the Up and Down Arrow buttons to highlight 10%, 20%, 50% or Continuous duty cycle.

When the buttons are pressed, an audible tone will be heard as the duty cycle toggles through the options on the LCD.

9. Press the Enter button to accept the selection.
10. Press the Down Arrow button to highlight Display.
11. Press the Enter button to select either Watts or W/cm². When the button is pressed, the display will toggle from Watts to W/cm² and back again as long as the button is being pressed.

12. Press the Up and Down Arrow buttons to select the manner in which you want to be notified when the sound head becomes uncoupled.

The following list details these options:

- Pause timer & beep - the timer stops and the unit beeps once.
- Pause timer, no beep - the timer stops but the unit does not give an audible tone.
- Run timer & beep - the timer continues to count down and the unit beeps once.
- Run timer, no beep - the timer continues to count down, but the unit does not give an audible tone.

NOTE: When the applicator becomes uncoupled during treatment, it continues to distribute ultrasound energy.

13. Press the Enter button to accept the selection.
14. Press the Time button and raise or lower treatment time using the up and down arrows.
15. Press the INTENSITY button and raise or lower the unit's output using the up and down arrows.
16. Press the Start button. The unit will beep 5 times and the ultrasound power will distribute the selected output.

NOTE: When treatment time has expired a tone will sound three times.

17. The therapy can be interrupted at any time by pressing the STOP or PAUSE buttons. When the STOP button is pressed, the applicator stops emitting ultrasound energy, and the unit returns to the default settings. To resume therapy, press the Start button.

When the PAUSE button is pressed, the timer pauses and the applicator stops emitting ultrasound energy. To resume therapy, press the PAUSE button again.

18. The parameters of the therapy can be changed at any time during the therapy session by pressing the appropriate button.
19. After therapy is complete, wipe excess ultrasound gel from the patient's skin and the sound head.

STARTING, STOPPING, AND INTERRUPTING COMBINATION THERAPY

Combination therapy consists of using electrotherapy and ultrasound therapy simultaneously. You may choose to use one or both channels of electrotherapy in conjunction with ultrasound.

Combination therapy utilizes the ultrasound modality in conjunction with High Voltage Pulsed Current (HVPC), VMS™, Symmetrical Biphasic and Asymmetrical Biphasic, Premodulated, or Interferential to generate a therapeutic effect. In this mode of therapy, the sound head of the ultrasound applicator becomes one half of the electrical circuit. An electrode attached to the black (-) Lead Wire completes the circuit.


To do this, do the following:

- Follow all appropriate procedures listed in the section entitled “**Electrotherapy Patient Preparation**” on page 26 and “**Ultrasound Patient Preparation**” on page 27.
- Connect the Lead Wires to the appropriate electrodes. For example, for all waveforms used by Channel 1 (all waveforms except Interferential), use the black (-) Lead Wire as the negative and the applicator as the positive. For the Interferential waveform, you will use three lead wires and the applicator for combination therapy: use the black (-) Lead Wire for Channel 1 as the negative electrode, the applicator as the positive, and the other Lead Wires for channel 2 as you would for an electrotherapy treatment.

NOTE: Do not use unnecessary force to connect the electrodes to the Lead Wires.
- Place the self adhesive electrodes on the sites prescribed by a qualified practitioner. Make sure you press them firmly on the patient’s skin to ensure good conductivity.
- Ensure the ultrasound applicator is plugged in.
- Depending on the type of waveform you intend to use, insert the Lead Wire into Channel 1 or both Lead Wire Connections on the Accessory Panel.
- On the Home screen, press the Up and Down Arrow buttons to highlight Combo.

- Press the Enter button. The Ultrasound parameters screen displays with the Edit Stim option highlighted.

NOTE: When you are using combination therapy for Channel 1 (all waveforms except Interferential), you may only set up electrotherapy for Channel 1. If you want to use both channels of electrotherapy in conjunction with ultrasound, you must first begin treatment. You may set up electrotherapy for Channel 2 later in these procedures (**step 21**).

- Press the Up and Down Arrow buttons to highlight Waveform.
- Press the Enter button. The Waveform menu displays.
- Press the Up and Down Arrow buttons to highlight the appropriate waveform.
- Press the Enter button to accept the selection. The Waveform menu closes.
- Press the Up and Down Arrow buttons to highlight Edit Stim.
- Press the Enter button to accept the selection. The current waveform’s parameters display.
- Press the Up and Down Arrow buttons to highlight the parameter you want to change and adjust it accordingly, making sure to press the Enter button after each adjustment.
- Press Enter to accept the parameters. The Parameter screen refreshes and the new parameters are displayed.
- Press the up and down arrows on the INTENSITY button to raise or lower the electrotherapy output.
- Press the Up and Down Arrow buttons to highlight Edit Ultrasound.
- Press the Enter button to accept the selection. The Ultrasound parameter screen displays.

STARTING, STOPPING, AND INTERRUPTING COMBINATION THERAPY (CONTINUED)


19. Continue with steps 4-23 of the section entitled **“Starting Stopping and Interrupting Ultrasound Therapy”** on pages 31-32.
20. Press the START button.
The unit will beep 5 times, you are returned to the Home screen, and the ultrasound and electrotherapy is distributed.
21. If you want to use Channel 2 for additional electrotherapy, continue with step 22. If you do not want to use additional electrotherapy, skip to step 25.
22. Press the Up and Down Arrow buttons to highlight Stim Channel 2.
23. Press the Enter button.
The Select Waveform screen displays.
24. Follow steps 8-12 of the section entitled **“Starting, Stopping, and Interrupting Electrotherapy”** on pages 29-30.
25. You may make parameter changes and stop or pause therapy on Channel 1, Channel 2, or Ultrasound by highlighting the appropriate form of therapy, and making the necessary changes.
26. When you have completed treatment, remove all electrodes from the patient.
27. Turn system power “OFF” by pressing the Power On/Off button. The unit beeps once and the blue light on the Power On/Off button flashes intermittently.

USING CLINICAL INDICATIONS


The indications contained in this section are to be used only as guidelines. Each patient should be individually assessed to determine the appropriateness of the parameter setting prior to use.

To select an indication for a patient, do the following:

1. On the home screen, press the Clinical Resources button. The Clinical Library window displays.
2. Using the Up and Down Arrow buttons, highlight Clinical Benefits and press the Enter button.

3. Using the Up and Down Arrow buttons, highlight the appropriate indication and press Enter. If a submenu displays, highlight the appropriate selection and press Enter.

4. Review the final parameters for the selected treatment. Make any necessary adjustments.
5. To begin therapy, prepare your patient according to the instructions outlined in the section entitled “Electrotherapy Patient Preparation” on Page 22 and Press Start.

CREATING A USER PROTOCOL FOR ELECTROTHERAPY


This is a library you create. You may store up to 15 protocols in the User Protocol Library. To create User Protocols for electrotherapy, do the following:

1. On the home screen, press the Up and Down Arrow buttons to highlight either Stim Channel 1 or Stim Channel 2.
2. Press the Enter button.


NOTE: User protocols can be used on any channel (Stim Channel 1, Stim Channel or Ultrasound). It does not matter on which channel they are created.

The Select Waveform screen displays.


3. Use the Up and Down Arrow buttons to highlight the appropriate waveform.

4. Press the Enter button. The Parameter screen displays.
5. On the parameter screen, press the Up and Down Arrow buttons to highlight the parameter you want to change and adjust it accordingly, making sure to press the Enter button after each adjustment.

6. Press the Clinical Library button. The Save Protocol screen displays.
7. Press the Up Arrow and Down Arrow buttons to highlight any unused user protocol. If you select the Save as Default protocol, this will become the protocol displayed when the waveform is selected on the Waveform screen.


8. Press the Enter button to accept the highlighted selection and save your custom protocol. The User Protocol Confirmation window displays to indicate that the protocol is now saved as the number you specified.
9. Press any button on the Operator Interface. The Parameter screen displays and your new protocol is now saved.

CREATING A USER PROTOCOL FOR ULTRASOUND THERAPY

To create User Protocols for ultrasound therapy, do the following:

1. On the Home screen, press the Up and Down Arrow buttons to highlight Ultrasound.
2. Press the Enter button.
The Ultrasound parameter screen displays.
3. On the parameter screen, press the Up and Down Arrow buttons to highlight the parameter you want to change and adjust it accordingly, making sure to press the Enter button after each adjustment.
4. Press the Clinical Library button.
The Save Protocol screen displays.
5. Press the Up Arrow and Down Arrow buttons to highlight any unused user protocol. If you select the Save as Default protocol, this will become the protocol displayed when Ultrasound is selected on the Home screen.

6. Press the Enter button to accept the highlighted selection and save your custom protocol.
The User Protocol Confirmation window displays to indicate that the protocol is now saved as the number you specified.
7. Press any button on the Operator Interface.
The Parameter screen displays and your new protocol is now saved.


CREATING A USER PROTOCOL FOR COMBINATION THERAPY

To create User Protocols for ultrasound therapy, do the following:

1. On the Home screen, press the Up and Down Arrow buttons to highlight Combo.
2. Press the Enter button.
The Ultrasound parameter screen displays.
3. Press the Up and Down Arrow buttons to highlight Waveform.
4. Press the Enter button.
The Waveform menu displays.
5. Press the Up and Down Arrow buttons to highlight the appropriate waveform.
6. Press the Enter button to accept the selection.
The Waveform menu closes.
7. Press the Up and Down Arrow buttons to highlight Edit Stim.
8. Press the Enter button to accept the selection.
The waveform's parameters display.
9. Press the Up and Down Arrow buttons to highlight the parameter you want to change and adjust it accordingly, making sure to press the Enter button after each adjustment.
10. Press the Up and Down Arrow buttons to highlight Edit Ultrasound.
11. Press the Enter button to accept the highlighted selection.
The Ultrasound parameter screen displays.
12. Press the Up and Down Arrow buttons to highlight the parameter you want to change and adjust it accordingly, making sure to press the Enter button after each adjustment.

CREATING A USER PROTOCOL FOR COMBINATION THERAPY (CONTINUED)

- Press the Clinical Library button.
The Save Protocol screen displays.


- Press the Up Arrow and Down Arrow buttons to highlight any unused user protocol. If you select the Save as Default protocol, this will become the protocol displayed when Combo is selected on the Home screen.
- Press the Enter button to accept the highlighted selection and save your custom protocol.
The User Protocol Confirmation window displays to indicate that the protocol is now saved as the number you specified.
- Press any button on the Operator Interface.
The Parameter screen displays and your new protocol is now saved.

USING PROTOCOLS

- On the home screen, press the Clinical Library button. The Clinical Library screen displays.
- Using the Up Arrow and Down Arrow buttons, highlight the Retrieve User Protocol option.
- Press the Enter button to accept the highlighted selection. A list of user-defined protocols displays.

- Press the Up and Down Arrow button to highlight the appropriate protocol. As you highlight each protocol, a description of the protocol's parameters displays to the right.

- Press the Enter button to select the highlighted protocol.
The Parameters screen displays the parameters of the protocol you selected.
- Verify the parameters of this program, and use the appropriate buttons on the Operator Interface to adjust any setting, if necessary. For example, to adjust the time, press the up and down arrows on the TIME button.
- To begin therapy, and depending on if you intend to apply electrotherapy, ultrasound therapy, or combination therapy, perform the appropriate procedures outlined on [pages 27-34](#).

SYSTEM UTILITIES

Changing Protocol Parameters

You may change any parameter prior to or during therapy. To make Intensity and Treatment Time changes, touch the respective buttons and use the up or down arrows to advance to the desired settings.

To change other parameters during therapy, do the following:


- On the home screen, use the Up and Down Arrow buttons to highlight the channel on which therapy is currently running.
- Press the Enter button to select the highlighted option.
The parameters of the current therapy session display.
- Press the Up and Down Arrow buttons to highlight the parameter you want to change and adjust it accordingly, making sure to press the Enter button after each adjustment.

SYSTEM UTILITIES (CONTINUED)

Changing Default Protocols for Electrotherapy

To change the power up presets, do the following:


1. On the home screen, press the Up and Down Arrow buttons to highlight either Stim Channel 1 or Stim Channel 2.
2. Press the Enter button to select the highlighted option.
3. Use the Up and Down Arrow buttons to highlight the appropriate waveform.
4. Press the Enter button.
The parameter screen displays.
5. Press the Up and Down Arrow buttons to highlight the parameter you want to change and adjust it accordingly, making sure to press the Enter button after each adjustment.
6. Press the Clinical Library button.
The Save Protocol screen displays.
7. Use the Up Arrow and Down Arrow buttons to highlight Save as Default protocol. This will become the protocol displayed when the waveform is selected on the Waveform screen.

8. Press the Enter button to accept the highlighted selection.
The Default Protocol Confirmation window displays.
9. Press any key to confirm the settings.
You are returned to the Clinical Library menu.

Changing Default Protocols for Ultrasound Therapy

1. On the home screen, press the Up and Down Arrow buttons to highlight Ultrasound.
2. Press the Enter button to select the highlighted option.
The Ultrasound parameters screen displays.
3. Press the Up and Down Arrow buttons to highlight the parameter you want to change and adjust it accordingly, making sure to press the Enter button after each adjustment.
4. Press the Clinical Library button.
The Save Protocol screen displays.


5. Use the Up Arrow and Down Arrow buttons to highlight Save as Default protocol. This will become the protocol displayed when Ultrasound is selected on the Home screen.
6. Press the Enter button to accept the highlighted selection.
The Default Protocol Confirmation window displays.
7. Press any key to confirm the settings.
You are returned to the Clinical Library menu.

Changing Default Protocols for Combination Therapy

1. On the home screen, press the Up and Down Arrow buttons to highlight Combo.
2. Press the Enter button to select the highlighted option. The Ultrasound parameters screen displays.
3. Press the Up and Down Arrow buttons to highlight Waveform.
4. Press the Enter button.
The Waveform menu displays.
5. Press the Up and Down Arrow buttons to highlight the appropriate waveform.
6. Press the Enter button to accept the selection.
The Waveform menu closes.

SYSTEM UTILITIES (CONTINUED)

7. Press the Up and Down Arrow buttons to highlight Edit Stim.
8. Press the Enter button to accept the selection.
The waveform's parameters display.
9. Press the Up and Down Arrow buttons to highlight the parameter you want to change and adjust it accordingly, making sure to press the Enter button after each adjustment.
10. Press the Up and Down Arrow buttons to highlight Edit Ultrasound.
11. Press the Enter button to accept the highlighted selection.
The Ultrasound parameter screen displays.
12. Press the Up and Down Arrow buttons to highlight the parameter you want to change and adjust it accordingly, making sure to press the Enter button after each adjustment.
13. Press the Clinical Library button.
The Save Protocol screen displays.
14. Use the Up Arrow and Down Arrow buttons to highlight Save as Default protocol. This will become the protocol displayed when Ultrasound is selected on the Home screen.


15. Press the Enter button to accept the highlighted selection.
The Default Protocol Confirmation window displays.
16. Press any key to confirm the settings.
You are returned to the Clinical Library menu.

Restoring Factory Protocols

If necessary, you can choose to restore the unit's original (default) waveform parameters when it was shipped to you.

NOTE: This procedure will erase all user-defined protocols. To restore the unit's original waveform parameters, do the following:

1. On the Home screen, press the Clinical Library button.
The Clinical Library window displays.
2. Press the Up Arrow or Down Arrow buttons to highlight Restore Factory Protocols.
3. Press the Enter button to accept the highlighted selection.
The unit displays the message "Restoring Protocols Please wait."
The user-defined protocols are erased and the waveforms are restored to the original parameters. Then the Restore Factory Protocols Confirmation window displays.


4. Press any button on the Operator Interface.
You are returned to the Clinical Library window.

SYSTEM UTILITIES (CONTINUED)


Changing Languages

To change the language displayed on the LCD, do the following:


1. Press the Clinical Library button.
The Clinical Library screen displays.

2. Use the Down Arrow and Up Arrow buttons to highlight the Language option.
3. Press the Enter button to accept the highlighted selection.
The Language menu displays.


4. Press the Up and Down Arrow buttons to highlight the appropriate language.
5. Press the Enter button to accept the highlighted selection.
The unit now displays the language you selected.

Restoring Factory Settings

To restore the original language on the unit, do the following:

1. On the Home screen, press the Clinical Library button. The Clinical Library screen displays.
2. Press the Up Arrow or Down Arrow buttons to highlight the Restore Factory Settings option.
3. Press the Enter button to accept the highlighted selection.
The Restore Factory Settings Confirmation screen displays.


4. Press any button on the Operator Interface.
The default power up settings are restored and you are returned to the Clinical Library screen.

Viewing Unit Version Information

Use this utility to determine the unit's software version. To do this, do the following:

1. On the Home screen, press the Clinical Library button. The Clinical Library screen displays.
2. Use the Up Arrow and Down Arrow buttons to highlight the View Unit Information option.
3. Press the Enter button to accept the highlighted selection.
The Unit Version Information screen displays.

4. Press any key to return to the Clinical Library window.

ERROR CODES

The Intellect Transport Combo displays error messages to inform the user of problems or potential problems with the unit, modality, or accessories. These are numbered so the user can possibly correct the problem without the aid of service personnel. Use the following Troubleshooting Charts to define the error codes, and locate the probable cause and possible remedies before contacting the dealer or factory for technical service.

Code Number	Type Message	Probable Cause	Possible Remedies
104	Message	User attempted to perform an electrotherapy session, but both channels are already in use.	A. Wait until the previous therapy session finishes. B. Press the STOP button to end the therapy session on either channel.
105	Message	User selected a two channel electrotherapy treatment, but at least one channel is already in use.	A. Wait until the previous therapy session finishes. B. Press the STOP button to end the therapy session on either channel.
107	Warning	Bad Contact Quality	A. Make certain Electrodes are making proper contact with the treatment area. B. Make certain Lead Wires are properly connected to Electrodes. C. Replace Electrodes and Lead Wires.
108	Warning	Shorted Lead Wires	A. Check Electrodes and Lead Wires. Make certain Lead Wires are not damaged and are properly connected to the system. Make certain Lead Wires are properly connected to the Electrodes and that electrodes are not damaged and are making proper contact with treatment area. B. Replace Lead Wires and Electrodes.
109	Message	While performing ultrasound therapy using the 10 cm ² sound at an intensity greater than 15 W (1.7 W/cm ²), user attempted to begin an electrotherapy session on a second channel.	A. Use a smaller sound head. B. Set the ultrasound at an intensity less than 15 W (1.7 W/cm ²). C. Wait until ultrasound session is complete.
200-399	Message	Refer to Warning on next page.	Refer to Warning on next page.

WARNING

In the event that an Error message or Warning appears beginning with a 2 or 3, immediately stop all use of the unit and contact the dealer or Chattanooga Group for service. Errors and Warnings in these categories indicate an internal problem with the unit that must be tested by Chattanooga Group or a Field Service Technician certified by Chattanooga Group before any further operation or use of the system. Use of a unit that indicates an Error or Warning in these categories may pose a risk of injury to the patient, user, or extensive internal damage to the system.

TROUBLESHOOTING ERROR MESSAGES

Troubleshooting Error Messages

The following messages are displayed on the LCD under the following conditions:

Message

Over Temp
Unplugged
Uncoupled

Displayed When

sound head reaches a temperature which could damage the crystal
sound head not plugged in or faulty sound head
sound head is not making good contact with the patient

REPLACEMENT AND OPTIONAL ACCESSORIES

The following provides users of the Intellect® Transport 2 the necessary information to order replacement accessories used with the system. This list of replacement accessories is designed for use with the Intellect® Transport 2. When ordering, provide the respective part number, description, and quantity desired.

STANDARD ACCESSORIES

Model Number	Description
47383	5 cm ² Ultrasound Applicator
47979	Transport 2 User Manual
48100	International Accessory Kit (for Combo device only)
40-0097	Leadwire Assembly, Ch 1 (for Combo device only)
40-0098	Leadwire Assembly, Ch 2 (for Combo device only)
79967-4	Silicon/Carbon Electrode 6x8cm Black (for Combo device only)
79970-4	Sponge Pocket 9x10.5cm (for Combo device only)
12-10648	Nylatex Wrap 2-1/2x24 Sewn (for Combo device only)
42197	Dura-Stick+ 5cm Round Electrodes (for Combo device only)

OPTIONAL ACCESSORIES

Model Number	Description
70010	Stim CH1/2 Leadwire Kit STD
10648	Nylatex Wrap (2 per pack)
79967-4	6 x 8 cm (2.5 x 3") carbon electrodes (4 per pack)
79967-40	6 x 8 cm (2.5 x 3") carbon electrodes (40 per pack)
79970-4	Sponge Pocket 9 x 10.5 cm (4 per pack)
79970-40	Sponge Pocket 9 x 10.5 cm (40 per pack)
47382	Transport 2 Ultrasound Applicator 2 cm ²
47383	5 cm ² Ultrasound Applicator
47384	Transport 2 Ultrasound Applicator 10 cm ²
14-1086	Battery
15-1136	Cart
14815	Main Power Cord UK
14816	Main Power Cord India
14819	Main Power Cord Australia
14820	Main Power Cord EU
15-3210	CARRYING BAG

CLEANING THE INTELECT® TRANSPORT 2

With the system disconnected from the power source, clean the system with a clean, lint-free soft cloth moistened with water and mild antibacterial soap. If a more thorough cleaning is needed, use a cloth moistened with an antimicrobial cleaner. Do not submerge the system in liquids. Should the unit accidentally become submerged, contact the dealer or DJO, LLC Service Department immediately.

Cleaning the electrotherapy lead wires

To clean the Lead Wires, disconnect them from the unit and wipe down with a clean, lint free cloth moistened with water and mild antibacterial soap .

Cleaning the Ultrasound applicator

To clean the ultrasound applicators, disconnect them from the unit and wipe down with a clean, lint free cloth moistened with water and mild antibacterial soap.

The ultrasound head may be disinfected with isopropyl alcohol between each therapy session.

MONTHLY INSPECTION

The following items should be checked at least monthly to ensure proper operation of this unit:

- Power cord and plug: Check to make sure the cord is not frayed, kinked, and does not have torn or cut insulation.
- Ultrasound applicator cable: Check to make sure the cable is flexible, free of kinks, not frayed, and the insulation is intact.
- Ultrasound head face: Check to see that there is no build-up of gel or foreign material on the aluminum face.

YEARLY CALIBRATION

The unit was calibrated during the manufacturing process and doesn't need calibration during the product life.

Annual factory calibration is required for all Ultrasound Applicators. Only the Applicators should be sent to the factory or a Trained Technician for this procedure.

DEVICE DISPOSAL

Council Directive 2012/19/EU concerning Waste Electrical and Electronic Equipment (WEEE) requires not to dispose of WEEE as municipal waste. Contact your local distributor for information regarding disposal of the unit and accessories.

ADDITIONAL INFORMATION

Go to the Chattanooga website www.chattanoogaarehab.com

WARRANTY REPAIR/OUT OF WARRANTY REPAIR

Service

When the Intellect® Transport 2 or any accessories require service, contact your selling dealer or your DJO Service Department contact.

Service to these units will be performed only by a service technician certified by the Company.

Expected Life

- Device expected life is five years
- Accessories expected life is one year
- Gel electrodes and ultrasound gel are shelf life accessories and their shelf life is less than device expected service life. Shelf life is indicated in electrodes packaging and on gel bottle.

WARRANTY

DJO LLC ("Company") warrants that the Intellect® Transport 2 devices ("Products") are free of defects in material and workmanship. This warranty shall remain in effect for two years (24 months) from the date of original consumer purchase. During the two-year warranty period from the date of delivery of the product to the end customer, defects will be remedied at no charge to the customer upon the customer furnishing adequate proof that the defect is due to defects in material or workmanship.

Attention

Modifications to the device are not permitted. Any unauthorized opening, repair or modification of the device by unauthorized personnel will relieve the manufacturer of its liability and responsibility for safe system operation. This will automatically void the warranty even before the end of the warranty period. The warranty period for accessories is 90 days. Accessories consist of Lead Wires and Electrodes. The warranty period for the Therapy System Cart and Ultrasound Applicators is one year (12 months).

This Warranty Does Not Cover:

- Replacement parts or labor furnished by anyone other than the Company, the selling dealer or a Company service technician
- Defects or damage caused by labor furnished by someone other than Company, the selling dealer or a Company service technician
- Any malfunction or failure in the Product caused by product misuse, including, but not limited to, the failure to provide reasonable and required maintenance or any use that is inconsistent with the Product User's Manual

COMPANY SHALL NOT BE LIABLE IN ANY EVENT FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES.

This warranty gives you specific legal rights and you may also have other rights which vary from location to location. The Company does not authorize any person or representative to create for it any other obligation or liability in connection with the sale of the Product.

Any representative or agreement not contained in the warranty shall be void and of no effect.

THE FOREGOING WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

ELECTROMAGNETIC COMPATIBILITY (EMC) TABLES

Device conforms to the following standards:	
AAMI STD ES 60601-1 CSA 22.2 NO. 6011 IEC 60601-1	Medical electrical equipment Part 1: General requirements for basic safety and essential performance
IEC 60601-1-6 NO. 60601-1-6	Medical electrical equipment – Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability
IEC 60601-2-5 NO. 60601-2-5	Medical electrical equipment - Part 2-5: Particular requirements for the basic safety and essential performance of ultrasonic physiotherapy equipment
IEC 60601-2-10 NO. 60601-2-10	Medical electrical equipment - Part 2-10: Particular requirements for the basic safety and essential performance of nerve and muscle stimulators

Guidance and manufacturer's declaration – electromagnetic emissions

The Intellect® Transport 2 is intended for use in the electromagnetic environment specified below.
The customer or the user of the Intellect® Transport 2 should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Intellect® Transport 2 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The Intellect Transport 2 is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded: Warnig: This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the Intellect Transport 2 or shielding the location.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

NOTE The EMISSIONS Characteristics Of this equipment make it Suitable for use in industrial areas and hospitals(CISPR 1 Class A). If it is used in a residential environment (for Which CISPR 11 Class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures. such as relocating or re-orienting the equipment.

ELECTROMAGNETIC COMPATIBILITY (EMC) TABLES (CONTINUED)

Guidance and manufacturer's declaration – electromagnetic immunity			
The Intelect® Transport 2 is intended for use in the electromagnetic environment specified below. The customer or the user of the Intelect® Transport 2 should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	<p>Risk assessment on the Intelect® Transport 2 indicates the compliance levels claimed are acceptable when ESD-precautionary measures are taken.</p> <p>The Intelect® Transport 2 may be susceptible to Electrostatic Discharge (ESD) at greater than ±7 kV when first grasping the Ultrasound applicator. In the event of such a discharge, the Intelect® Transport 2 may display a permanent error. The Intelect® Transport 2 will terminate all active outputs (stim, ultrasound), automatically place the unit in a safe state.</p> <p>To prevent Electrostatic Discharge (ESD) at greater than ±7 kV:</p> <ul style="list-style-type: none"> • Grasp and hold the Ultrasound applicator prior to starting treatment. If the applicator must be put down prior to completion of treatment, stop the current treatment first and then place the applicator in the holder. • Maintain humidity in the use environment to at least 50% relative humidity. • Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, DJO recommends implementing additional controls to maintain relative humidity to at least 50%. • Communicate these ESD-precautionary procedures to healthcare staff, contractors, visitors and patients.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) +/- 2 kV line(s) to ground	± 1 kV differential mode +/- 2 kV line(s) to ground	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles Single phase at 0°	0% UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles Single phase at 0°	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Intelect® Transport 2 requires continued operation during power mains interruptions, it is recommended that the Intelect® Transport 2 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/M	30 A/M	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Proximity magnetic fields IEC 61000-4-39	CW; 8 AM 2.1 KHZ; 65 A/M 50 KHZ; 7.5 A/M	CW; 8 AM 2.1 KHZ; 65 A/M 50 KHZ; 7.5 A/M	Test frequency 30 kHz Test frequency 134.2 kHz Test frequency 13.56 MHz
NOTE UT is the a.c. mains voltage prior to application of the test level.			

ELECTROMAGNETIC COMPATIBILITY (EMC) TABLES (CONTINUED)

Guidance and manufacturer’s declaration – electromagnetic immunity			
<p>The Intellect® Transport 2 is intended for use in the electromagnetic environment specified below. The customer or the user of the Intellect® Transport 2 should assure that it is used in such an electromagnetic environment.</p>			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz outside ISM bands ^a	3 Vrms	The device is intended for use in a typical professional healthcare facility environment. Portable RF communications equipment should be used no closer than 30 cm (12 inches) to any part of the Intellect® Transport 2, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
	6 Vrms 150 kHz to 80 MHz in ISM bands ^a	6 Vrms	
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	
	9-28V/m in wireless bands	9-28V/m	
			
<p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p>			
<p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p>a) The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz.</p>			

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DJO, LLC, 5919 Sea Otter Place, Suite 200, Carlsbad, CA 92010

[enovis.com/uk](https://www.enovis.com/uk)

Individual results may vary. Neither Enovis, DJO, LLC or any of their subsidiaries dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.

11/2024
PN 47979 – Rev D

Copyright © 2024 DJO, LLC, a subsidiary of Enovis Corporation

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Copyright © 2024 DJO, LLC, subsidiaria de Enovis Corporation
Copyright © 2024 DJO, LLC, une filiale d'Enovis Corporation
DJO France SAS
MDSS CH GmbH